


HSW ECO-MATIC[®] TWIN

The Market's First Light Weight, Ergonomic Compact Double Barrel Plastic Syringe

- Simultaneous vaccination of two products independently
- Premium-quality plastic syringe with durable plastic barrels
- Reusable-The only plastic syringe that can be boiled
- Ergonomic, lightweight design for comfort and ease of use
- Bottle mount and/or tubing connection
- Adjustable barrel orientation for comfort and different bottle sizes



- Unique Patent Pending Needle Nut Swivel Design

Needle distance can be adapted to the age, size and species of the animal.
Seven different distances:

20 mm – 21.1 mm – 24 mm – 28 mm – 32 mm – 34.9 mm – 36 mm

- Four Fixed Dose Barrel combinations available:

Package Includes:

- Two Universal bottle collars (20 mm & 30 mm)
- Two Universal spigot draw off collars (4 - 6 mm)
- Two soft PVC tubing

Allflex SKU	Allflex Description
1/1 EM-TWIN	1/1 EM TWIN Syringe With Accessory Pack
2/2 EM-TWIN	2/2 EM TWIN Syringe With Accessory Pack
1/2 EM-TWIN	1/2 EM TWIN Syringe With Accessory Pack
.5/2 EM-TWIN	.5/2 EM TWIN with Accessory Pack

800.989.TAGS • allflex.global/us

Allflex Livestock Intelligence • P.O. Box 612266 • Dallas/Fort Worth Airport, TX 75261-2266

ALIDET20_001