



- 2ml and 5ml sizes available
- Automatic self filling with bottle attachment
- Heavy-Duty Bottle Cage for bottle protection (50ml & 100ml bottles only)
- 2EM-BT adjustable dosage from .25ml to 2ml
- ♦ 5EM-BT adjustable dosage from .5ml to 5ml
- Metal Luer Lock Tip
- ♦ Removable barrel for cleaning
- Durable Plastic Design
- Easy V-grip for high volume use
- ◆ Amber Colored Barrel reduce UV sunlight exposure "Allflex Exclusive Feature"





P.O. Box 612266 * Dallas/Fort Worth Airport,TX 75261-2266 800-989-TAGS (8247) * Customer Toll Free FAX 877-456-3969



Allflex Canada

4135 Berard * St. Hyacinthe, Quebec J2S 8Z8 866-505-TAGS (8247) * FAX 450-261-8028