



# Allflex® Tissue Sampling Unit (TSU)

PRESERVING A WEALTH OF GENETIC INSIGHTS



## ALLFLEX TSUS ARE YOUR CONVENIENT GATEWAY INTO GENETIC IMPROVEMENT

- One quick squeeze of the Allflex TSU Applicator deposits a clean, consistent DNA sample into a liquid preservative.
- Inside that vial resides a wealth of genetic insight on that animal's pedigree and impact on future progeny.
- Combined with your expertise, DNA testing can speed progress on goals and fill your herd with superior cattle.

## BETTER INSIGHTS DRIVE BETTER HERD MANAGEMENT DECISIONS

DNA reports help you firm up the decisions that influence success in herd improvement. Is a calf a true progeny of a sire? Are all the bulls doing their job? Which heifers should be kept in the breeding herd? What are the best sire genetics for your herd? How do cows score on indexes?

## FASTER, CLEANER AND EASIER - A GOOD TEST BEGINS WITH A GOOD DNA SAMPLE

Allflex TSUs are faster, easier and more effective for you. They also are the industry's preferred DNA sampling method, preventing cross contamination or messy sampling errors common with blood or hair.





# Streamlining your use of production records

## MAKING TISSUE SAMPLING MORE CONVENIENT

With TSUs you can take samples when processing cattle and decide later about lab testing. For record keeping, each sterile vial is individually numbered and barcoded for storing samples, ordering tests or comparing DNA test results to other animal information.



The TSU Applicator deposits a tissue sample into the sterile vial for fast, clean, easy DNA sampling.



## MATCHING TSU VIALS WITH OTHER FORMS OF ANIMAL ID

You can purchase matched sets linked by the same identifier number on a TSU, EID or visual tag. Matched sets make it easier to link animal ID numbers to production or progeny records.

Solutions Available for Tissue Collection and DNA Testing	
<p><b>Allflex TSU Applicator</b> Applicator for collecting TSU samples SKU: APP-TSU • UIN: 321125</p>	
<p><b>TSU 10-Pack</b> Ten barcoded TSUs pre-numbered or custom numbered. SKU: AF-TSU-10 • UIN: 321123</p>	
<p><b>TSU Matched Sets with Visual and/or EID Tag</b> Multiple combination options available. Shown at right as matched tamperproof tag, EID and TSU. SKU: TSU-HDX840/GTLF-Y-SET1 • UIN: 320054</p>	
<p><b>100-pack TSUs, Pre-numbered or Blank</b> 100 vials, pre-numbered or custom numbered. SKU: AF-TSU-100 (pre-numbered) • UIN: 321124</p>	
<p><b>Allflex TSU Shipping Trays</b> Tray holds up to 96 vials for protected storage or shipping. SKU: TSU-TRAY+LID • UIN: 310399</p>	

## Tagging Today. Tagging Tomorrow.

[www.allflexusa.com](http://www.allflexusa.com)

Copyright ©2025 Merck & Co., Inc., Rahway, NJ, USA and its affiliates. All rights reserved.  
US-Allflex-251200007

