



COX & RITCHEY TRADE PARTNER BROCHURE

2017/18





Welcome to the 2017/18 Cox & Ritchey Trade Partner Brochure which showcases our wide range of livestock consumables, stock management equipment and livestock identification tags. As exclusive UK distributors for the Swiss made Heiniger clipping and shearing equipment you'll find some of the range, which is widely acknowledged for its precision engineering, build quality and reliability on pages 50-56.

Comprehensive marketing support for your business

Cox & Ritchey are committed to supporting you throughout the course of the year whether this is with point of sale marketing material, specific product training or product support literature.

Our Account Managers are also happy to provide practical advice and demonstrations for all our **Prattley** and **Alligator** Sheep and Cattle Handling Systems and our **Tru-Test** Weighing Equipment. This of course is not an exhaustive list so if there is something which you would like our help with to promote please do not hesitate to contact us.

Dedicated Customer After Sales Service

Finally, to complement some of the products contained within the brochure we also provide a dedicated Customer After Sales Service which includes an in-house servicing and repair facility for **Tru-Test** Electronic Weighing Equipment and EiD Readers and **Heiniger**, **Oster** and **Andis** branded clipping and shearing equipment.

We also offer in-house servicing on selected Gas and Electric Dehorners (as sold by us), **Hotshot** Pulsar Coaxers and **Philips** Drenchers and Vaccinators as well as a comprehensive blade sharpening service all of which are an essential part of the customer aftercare service provided by **Cox & Ritchey**. More details of our Customer after Sales Service can be found on page 58.

The Cox & Ritchey team



CONTENTS

Introduction	02
Birthing	04
Nutrition	06
Calf Feeders	08
Calf and Cow Care	10
Cattle Breeding	13
Drenching and Vaccinating	15
Foot Care	20
Lambing	22
Mating and Marking	35
Sundries	38
Animal Handling and Weighing Equipment	40
Heiniger Clipping and Shearing	50
Grooming	57
Customer After Sales Service	58
Livestock Identification	59
Index	70

HOW TO ORDER

You can place your order in a number of ways either by contacting our customer services team on:

 **0845 6008081**

 **tradesales@allflex.co.uk**

 **0800 7836655**

Or alternatively please contact your Cox & Ritchey Account Manager:

Your Account Manager

Mobile No.

Email Address

Your Account No.

POSTAGE CHARGES

Tags

Part Code	Product Type	Description	Cost
POST1	Tags	Replacements	£1
POST2	Tags	All other Tags	£2

Products

Part Code	Product Type	Description	Cost
POST3	Stock	Under £150 (excl. those with special pricing)	£6
POST4	Stock	Over £150	FOC

Special Deliveries

Part Code	Product Type	Description	Cost
POST7	Special	Royal Mail next working day guaranteed by 9am	£15
POST8	Special	Royal Mail next working day guaranteed by 1pm	£7
POST9	Special	UK Mail next working day delivery pre-10.30am*	£17.18
POST10	Special	UK Mail next working day delivery pre 12pm	£11.27
POST11	Special	UK Mail Saturday delivery (any time) *	£10
POST13	Special	UK Mail Saturday delivery pre 9am*	£22.82
POST14	Special	UK Mail Saturday delivery pre-10.30am*	£20

PLEASE NOTE: Postal charges are subject to change without notice.


RITCHEY THICK SUPERLUBE



Ritchey Superlube is a thick general purpose gel used for animal husbandry it contains a preservative which is non-spermicidal. The gel is ideal for use in the delivery of new born animals.

PART NO	DESCRIPTION
52801	Superlube Thick 500 ml - Pack Size 12
52802	Superlube Thick 1L - Pack Size 12
52803	Superlube Thick 2.5L - Pack Size 4
52804	Superlube Thick 5L - Pack Size 2
52805	Superlube Thick 25L - Pack Size 1

HK HEAVY DUTY CALVING AID



The HK Heavy Duty Calving Aid is an essential tool for difficult calving and can help reduce the number of calves lost at birth.

The calving aid is supplied with a galvanised steel two piece 1800mm shaft and a pair of ropes. The aid comes with a 20-20 ratchet which allows for alternate or straight pulls and includes a simple easy-to-use quick-release mechanism even under high loading.

KEY FEATURES INCLUDE:

- Lightweight and easy to use single person operation
- Simple ratchet action
- Comes complete with two nylon calving ropes
- The buttock is made with hard moulded plastic

PART NO	DESCRIPTION
54012	HK Heavy Duty Calving Aid
54013	Pair of Spare Ropes for HK Calf Puller

* Spares available

CALVING ROPE



An essential aid used during calving and manufactured using rot resistant nylon for prolonged use.

PART NO	DESCRIPTION
433019	Calving Rope 1 Loop
433020	Calving Rope 2 Loops

NITRILE GLOVES

A synthetic rubber glove which provides tough chemical resistance coupled with maximum sensitivity. Available sizes: Small, Medium, Large and Extra Large.

KEY FEATURES INCLUDE:

- Provides tough chemical resistance with maximum sensitivity
- Reduces the risk of contamination
- Latex free, non-powdered, non-sterile, contains no protein

PART NO	DESCRIPTION	PART NO	DESCRIPTION
496305S	Nitrile Gloves Pk 100 Small	496305L	Nitrile Gloves Pk 100 Large
496305M	Nitrile Gloves Pk 100 Medium	496305XL	Nitrile Gloves Pk 100 Extra Large



DISPOSABLE GLOVES

These high quality one size disposable gloves are ideal for food preparation and general purpose use.

The gloves are smooth which allows for easier application and come packaged in an easy to use dispensing box. These gloves can also be used for catering, cleaning and can be used by mechanics.

- Available in packs of 100

PART NO	DESCRIPTION
54021	Disposable Gloves 100pk - Not Illustrated

SILKTOUCH ARM LENGTH GLOVES

Strong and robust arm length gloves, which provide maximum sensitivity making these gloves ideal for allergy sufferers.

- Available in packs of 100

PART NO	DESCRIPTION
496261	Silktouch Arm Length Gloves



ARM LENGTH DISPOSABLE GLOVES

Ideal for lambing, cleansing, AI, or where a high degree of hygiene is required. These sensitive watertight gloves fit either hand.

- Available in packs of 100

PART NO	DESCRIPTION
496133	Arm Length Disposable Gloves



LATEX GLOVES

KEY FEATURES INCLUDE:

- Strong and robust
- Ideal for short term use / frequent changing
- Available in sizes: Small, Medium, Large and Extra Large

PART NO	DESCRIPTION	PART NO	DESCRIPTION
496297S	Latex Gloves Box 100 Small	496297L	Latex Gloves Box 100 Large
496297M	Latex Gloves Box 100 Medium	496297XL	Latex Gloves Box 100 Extra Large



FARMER'S CHOICE CALF COLOSTRUM*

Farmers Choice Calf Colostrum is an EBL free, easy to mix, highly stable complementary feed for new born calves that are deprived of or are deficient in colostrum from their mothers at birth. The colostrum should be fed to all calves requiring colostrum and/or extra energy, preferably in the first 6 hours of life.

Weak calves and those where the mother's colostrum is suspect will benefit for up to 48 hours after birth.

FARMER'S CHOICE



KEY FEATURES INCLUDE:

- Actively supports calves during their critical new born stage
- Scientifically formulated to provide instant energy whilst supporting the calf's active immunity
- The natural nucleotides offer additional support to the immune system and gut function
- Simply add water and shake to mix
- Stays mixed so it is easy to feed, even in the harshest conditions
- Convenient single feed container

PERFECT TO USE WITH THE SPEEDY DRENCHER
See page 9

PART NO	DESCRIPTION
720 001	Farmer's Choice Calf Colostrum (Single Feed - 150g)
720 005	Farmer's Choice Calf Colostrum (5 Feed Pail - 750g)

* Sold as a single item

Also available...
LAMB COLOSTRUM
See page 25

FARMER'S CHOICE OVINE CONDITIONER



Farmer's Choice Ovine Conditioner provides a quick and effective way to deliver high levels of essential chelated trace elements, vitamins and amino acids to improve the condition of all sheep.

The conditioner incorporates organically chelated trace elements which stay in the body longer, providing a ready supply of trace elements, as and when the animal requires them.

KEY FEATURES INCLUDE:

- Quick and effective
- Incorporates organically chelated trace elements which stay in the body longer

PART NO	DESCRIPTION
710250CU	Farmer's Choice Ovine Cond with Copper 5Lt - Pack of 2
710225	Farmer's Choice Ovine Cond without Copper 2.5Lt - Pack of 4

DOSAGE RATES

EWES AND RAMS PRE-TUPPING

Breed	Measure	Timings
Rams	20ml	6-8 weeks before turnout & again at turnout
Hill ewes & small breeds	15ml	3-4 weeks prior to ram turnout
Lowland ewes & larger breeds	20ml	3-4 weeks prior to ram turnout
Intensively flushed lowland ewes	15ml	3 weeks prior to ram & 15ml at ram introduction

PRE-LAMBING

Breed	Measure	Timings
Hill ewes & small breeds	15ml	4-6 weeks pre-lambing
Lowland ewes & larger breeds	20ml	4-6 weeks pre-lambing

LAMBS & HOGGS

Breed	Measure	Timings
Young lambs	5ml	At marking or first worming
Lambs at 3 months or at weaning	10ml	Repeat every 3-4 months
Hoggs & ram lambs	10ml	Repeat every 3-4 months

FARMER'S CHOICE COBALT BOLUS

Farmer's Choice Cobalt Boluses are a long acting, dietetic complementary feed. Powerful boluses designed to provide sheep with their daily cobalt requirements for a full year.

Annual dosing is more convenient, requires less animal handling and gives peace of mind, knowing the flock's cobalt requirements are met.

Farmer's Choice boluses are small and powerful, measuring just 2cm long, and weighted to remain in the reticulum. They are easy to administer to both growing lambs and breeding stock.

A simple to use, quality constructed bolus gun is available to administer the boluses.

KEY FEATURES INCLUDE:

- Supplied in packs of 100 boluses
- Feed: 1 bolus, per animal, per year
- Release rate: 1mg Co per day



PART NO	DESCRIPTION
446014	Farmer's Choice Cobalt Bolus Pack Of 100

BALLING GUN



Sheep Balling Gun for administering Farmer's Choice Cobalt boluses. Specially designed to protect throat.

PART NO	DESCRIPTION
445071	Sheep Balling Gun

SHEEP RUBBER MOUTH PIECE

PART NO	DESCRIPTION
445083	Sheep Rubber Mouth Piece - Not illustrated

SUPERSTART CALF FEEDING SYSTEM



The Superstart calf feeding bottle with teat is a robust traditionally shaped feeding bottle which can be fitted with a stomach feeding tube for administering fluids to weakly calves.

PART NO	DESCRIPTION
504402	Superstart Bottle and Teat - 2.5 litres capacity
504404	Superstart Teat – Pack of four
504406	Superstart Cap
504408	Superstart Feeding Tube

SUPERSTART CALF TEAT



Superstart rubber teats (for use with the Superstart Calf Feeding System) have similar texture and 'feel' to calf bar feeders, to encourage easier transfer from the bottle to the bar.

KEY FEATURES INCLUDE:

- Encourages easier transfer from bottle to the bar

PART NO	DESCRIPTION
504404	Superstart Calf Teat (Pack of 4)

SUPERSTART STOMACH TUBE



KEY FEATURES INCLUDE:

- Easily fitted to the superstart feeding bottle
- Allows for easy administration of fluids to weakly calves

PART NO	DESCRIPTION
504408	Superstart Stomach Tube

CALF TEAT

A rubber calf teat which can be pulled over a bottle neck for feeding.



PART NO	DESCRIPTION
505067	Calf Teat

CALF TEATS



A supple, yet extremely strong pair of rubber white soft and medium red teats.

PART NO	DESCRIPTION
56201	Calf Teat - Red - Packs of 5
56202	Calf Teat - White - Packs of 5

CALF BUCKET & TEAT



A sturdy calf bucket with teat and handle which has a flat back to allow it to hang securely against hurdles.

PART NO	DESCRIPTION
504315	Calf Bucket and Teat
504129	Spare Teat

5L BLUE FEEDING BUCKET

KEY FEATURES INCLUDE:

- Flexible and robust
- 1L measuring increments
- Stainless steel handle with rotational grip



PART NO	DESCRIPTION
56180	5L Blue Feeding Bucket

FLEXIBLE ORAL CALF FEEDER*



The flexible oral calf feeder aids against the onset of dehydration and is made from non-toxic plastic. The oesophageal tube is moulded as one unit and designed to not enter the trachea.

PART NO	DESCRIPTION
504208	Flexible Oral Calf Feeder - 2.5 litre capacity

* Recommended single use only

RITCHEY EMERGENCY FEEDER



This emergency feeder is a must have item to have around the farm. Administer emergency colostrum or glucose directly into the stomach of newly born calves. Comes complete with a zip up bag and stomach tube.

KEY FEATURES INCLUDE:

- Zip up bag makes it easier to use
- Stomach tube is included

PART NO	DESCRIPTION
56000	Ritchey Emergency Feeder - Calf 2.5 litre capacity

SPEEDY FEEDER

A three speed milk feeder designed to make feeding quick and easy.

KEY FEATURES INCLUDE:

- A Low off centre teat for easier feeding and full draining
- 3 Speed pump action system ensures that weak or premature calves receive a constant flow of milk
- Ergonomically designed for front and back facing feeding

PART NO	DESCRIPTION
56231	Speedy Feeder - 2.5 litre capacity
56294	Speedy Feeder Teat and Top



SPEEDY DRENCHER

The Speedy Drencher is ideal for administering colostrum to new born calves.

KEY FEATURES INCLUDE:

- Ergonomic bi-directional handle bar which provides a better grip
- Practical screw top with large thread
- Wide base for maximum stability
- Drencher probe length 48cm with a capacity to hold 2.5L

PART NO	DESCRIPTION
56232	Speedy Drencher Complete
56295	Speedy Drencher Spare Probe




Perfect to use with Farmer's Choice Calf Colostrum. See page 6

LEVER ACTION DEHORNING SHEAR

For young cattle.

KEY FEATURES INCLUDE:

- Blade opening 57mm (2.25")
- Concave blades for a closer cut. 58cm (23")



PART NO	DESCRIPTION
410226	Lever Action Dehorning Shear

GEARED DEHORNING SHEAR

For young cattle.

KEY FEATURES INCLUDE:

- Powerful geared action.
- Blade opening 38mm (1.5")
- Concave blades for a closer cut 58 cm (20")



PART NO	DESCRIPTION
410093	Geared Dehorning Shear

SPARE BLADE SET

Spare blades for Lever Action Dehorning Shears
Part Code 410226

PART NO	DESCRIPTION
410457	Spare Blade Set - Not Illustrated

DE-BUDDER WITHOUT REGULATOR

A lightweight and easy to use Gas De-budder.

KEY FEATURES INCLUDE:

- Fully portable it can be used anywhere
- Consistent and constant heat capable of dehorning 20-25 calves in just one hour
- Quick – Averaging a speed of 10-12 seconds per horn
- Simple design, easy to maintain and easily repairable
- Economical to run, uses either Butane or Propane gas
- Reaches its working temperature in just 2-3 minutes
- Low maintenance, with screw fit brass head and grub screw



PART NO	DESCRIPTION
410135	De-budder without Regulator
410172	Propane Regulator
410170	Butane Regulator

PULSAR ELECTRONIC LIVESTOCK PROD

An electronic Livestock prod designed to provide assistance for the movement of cattle.

PART NO	DESCRIPTION
70020	Pulsar Fx22 Flexible Shaft - shaft size 12" - Not Illustrated
70021	Pulsar Fx32 Flexible Shaft - shaft size 22"
70022	Pulsar Fx42 Flexible Shaft - shaft size 32"

COPPER DEBUDDER FIRE TYPE

The Fire Type comes with a large copper head for rapid heating, no trailing wires or tubing. It is extremely easy to handle, which means less stress on the animal. No maintenance is needed and is ideal where there is no access to mains electricity .



KEY FEATURES INCLUDE:

- Easy to handle
- No maintenance required
- No trailing wires or tubing


PART NO	DESCRIPTION
409078	Copper De-budder Fire Type

BARNES TYPE DEHORNER (LARGE JAW)

Suitable for more mature cattle.

KEY FEATURES INCLUDE:

- Lightweight (1.44kg) for easy handling
- Blade opening 70mm



PART NO	DESCRIPTION
410202	Barnes Type Dehorner (Large Jaw)

BARNES TYPE DEHORNER

A lightweight and easy to use Dehorner suitable for use on calves aged between 2-6 months.

KEY FEATURES INCLUDE:

- Blade opening 60mm. 33cm (13")
- Lightweight and easy to use



PART NO	DESCRIPTION
410111	Barnes Dehorner

CALF DEHORNER MK III PORTASOL



A fully portable and extra durable dehorner suitable for dehorning young calves between 2-8 weeks old. The Dehorner reaches its working temperature rapidly and features both click to ignite and manual ignition, ensuring fast, clean and low stress dehorning.

KEY FEATURES INCLUDE:

- Copper tip for rapid heat up to 650°C / 1200 F
- Extra wide tip available for dehorning older calves
- Click to ignite plus manual ignition
- Comes with 12 month warranty
- Also available: Gas Tank & Burner Assembly - MKIII

PART NO	DESCRIPTION
56503	Calf Dehorner MkIII Portasol
56508	Standard Spare Tip - 15mm
56507	Large Spare Tip 18.5mm
56514	Gas Tank & Burner Assembly - MKIII

HEAVY DUTY DEHORNING SHEARS



For mature cattle these dehorners have a powerful guillotine action and are made with high quality steel blades.

Blade opening 89 mm (3.5"). Length 1.2m (4') 10.5kg.

PART NO	DESCRIPTION
410068	Dehorning Shears

PLEASE NOTE: This is a large heavy high quality instrument, please check weight and dimensions before placing your order.

SPARE BLADE SET

Spare blades for 410068 Heavy Duty Dehorning Shears.

PART NO	DESCRIPTION
410081	Spare Blade Set - Not Illustrated

LEISS DEHORNING WIRE



The Leiss Dehorning Wire is a quick and easy to use steel wire used for dehorning adult cattle.

KEY FEATURES INCLUDE:

- Flexible and robust
- 12m / 37'

PART NO	DESCRIPTION
410056	Leiss Dehorning Wire

LIGHT GAUGE DEHORNING WIRE

Also used for embryotomy 10.5m (35').

PART NO	DESCRIPTION
410020	Light Gauge Dehorning Wire - Not Illustrated

HEAVY GAUGE DEHORNING WIRE

10.5m (35').

PART NO	DESCRIPTION
410032	Heavy Gauge Dehorning Wire - Not Illustrated

PAIR DEHORNING WIRE HANDLES



A pair of wire handles to safely use with dehorning wire.

PART NO	DESCRIPTION
410044	Pair Dehorning Wire Handles

STANDARD CALF JACKET



The standard calf jacket provides a water-resistant layer for new calves helping to avoid cold stress. The jacket is easy to fit, using simple nylon straps and clip fasteners. It provides comfort and protection to a farmer's valuable livestock.

KEY FEATURES INCLUDE:

- Stitched with rot resistant Kevlar thread for a longer lifespan
- 10cm wide / 58cm long
- Washable

PART NO	DESCRIPTION
56806	Standard Calf Jacket

DOWNACOW HARNESS



The Downacow harness maximises the recovery rate of a downed cow, whilst requiring the minimum amount of farm labour. Fabricated from rot resistant, nylon webbing. The Downacow Harness can be easily and quickly fitted by two people, and fully supports the cow, whilst remedial action is taken. The Downacow Harness is an economic solution to a troublesome problem.

KEY FEATURES INCLUDE:

- Maximises the recovery rate of the cow
- Constructed from rot resistant webbing

Maximum weight 500kg

PART NO	DESCRIPTION
404226	Downacow Harness

COX HEAVY DUTY HOBBLER



An extremely substantial set of cattle hobbles made from broad, soft feel nylon cuffs, which avoids abrasion on the heels of the animal. The hobbles are linked by a strong chain, incorporating a swivel. The hobbles are stitched with rot resistant Kevlar thread for a longer lifespan on the product.

KEY FEATURES INCLUDE:

- Stitched with rot resistant Kevlar thread
- Longer lifespan
- 10cm wide / 58cm long

PART NO	DESCRIPTION
404202	Cox Heavy Duty Hobbles

COPPER BULL RING



The Bull Ring is made from high quality copper and has a self-piercing point and comes complete with a brass shear screw.

KEY FEATURES INCLUDE:

- Made from high quality copper
- Self-piercing point

PART NO	DESCRIPTION
401031	Copper Bull Ring 64mm (2.5")
401043	Copper Bull Ring 70mm (2.75")
401055	Copper Bull Ring 75mm (3")
401122	Copper Bull Ring 100mm (4")

BULL RING APPLICATOR



For use with 75mm (3") rings.

PART NO	DESCRIPTION
401171	Bull Ring Applicator

HEAVY DUTY COW IMMOBILIZER

Fully adjustable and easy to fit, prevents nervous animals from kicking.

PART NO	DESCRIPTION
404305	Heavy Duty Cow Immobilizer - Not Illustrated

COW IMMOBILIZER TIP

PART NO	DESCRIPTION
404408	Cow Immobilizer Tip - Not Illustrated

COW IMMOBILIZER SPRING

PART NO	DESCRIPTION
404425	Cow Immobilizer Spring - Not Illustrated

MILK SUCKING PREVENTORS

WILSON'S TYPE MILK SUCKING PREVENTOR



Polished aluminium.

PART NO	DESCRIPTION
502054	Wilson's Type Milk Sucking Preventor

SPIKED METAL MILK SUCKING PREVENTOR



PART NO	DESCRIPTION
502078	Spiked Metal Milk Sucking Preventor

SPIKED PLASTIC MILK SUCKING PREVENTOR



PART NO	DESCRIPTION
502108	Spiked Plastic Milk Sucking Preventor

HALTERS

COTTON CALF HALTER



Adjustable white rope.

- Length 2.42m
- Diameter 1.1cm
- Circumference 2.54cm

PART NO	DESCRIPTION
406119	Cotton Calf Halter

COTTON CATTLE HALTER



Adjustable with metal ring, white rope.

- Length 3.1m
- Diameter 1.4cm
- Circumference 3.8cm

PART NO	DESCRIPTION
406090	Cotton Cattle Halter

CATTLE HALTER

White cotton webbing and rope.

- Length 2.51m
- Diameter 1.2cm
- Circumference 3.8cm

PART NO	DESCRIPTION
406077	Cattle Halter - Not illustrated

CATTLE HALTER, NYLON ROPE

Adjustable halter with metal ring, white rope.

- Length 2.81m
- Diameter 1.3cm
- Circumference 3.8cm

PART NO	DESCRIPTION
406107	Cattle Halter, Nylon Rope - Not illustrated

NYLON SHEEP / CALF HALTER

Adjustable white rope.

- Length 2.45m
- Diameter 1.2cm
- Circumference 1.6cm

PART NO	DESCRIPTION
406120	Nylon Sheep / Calf Halter - Not illustrated



PHILLIPS DRENCHING PRODUCTS

The NJ Phillips range of drenching instruments, are manufactured to the highest quality using the most up-to-date range of materials. This ensures complete compatibility with modern animal health formulations. With a heritage which dates back to 1930, the Phillips brand is based on quality and accuracy in manufacture, working closely with the agricultural industry and livestock producers they have built up a range of drenchers renowned for outstanding industry dosing accuracy.

2ML STF BOTTLE MOUNT INJECTOR



The bottle Mount Injector can be attached directly to the syringe, there is no need for tubes with this product. There is 0.2ml graduations from 0.2 - 2ml. The plastic Luer tip fits most Luer specification needles.

KEY FEATURES INCLUDE:

- Screw thread adjuster.
- Ergonomic 'finger-grip' design

PART NO	DESCRIPTION
40501	2ml Stf Bottle Mount Injector

PHILLIPS REPEATER SYRINGE



The Phillips Repeater Syringe is easy to use and hassle free. Fill direct from the bottle, set the dose and start injecting. When injecting ewes and lambs the dose is altered easily with a simple twist of the dose adjuster. No trailing tubes or awkward bottle mounts. Dial-a-dose setting can be altered at any time. Accepts disposable and reusable luer lock needles.

KEY FEATURES INCLUDE:

- Works well with all consistencies of medication, thick, thin and oily
- Spare needle holder/magazine in the handle
- Supplied with needle spanner
- Particularly useful for larger doses from small bottles

Available in two different sizes 25ml and 50ml

PART NO	DESCRIPTION
529175	Phillips Repeater Syringe 25ml
529606	Phillips Repeater Syringe 50ml

PHILLIPS SEMI. DISP. VACCINATOR 5ML WITH DRAW OFF



This high quality precision Phillips instrument is one of the most popular automatic plastic vaccinators in the UK. Lightweight, and economic, giving accurate administration of vaccine. The vaccinator comes complete with a draw off system and the injector is used as part of the popular Sterivax system.

KEY FEATURES INCLUDE:

- High quality precision
- Lightweight and economic
- Allows accurate administration of the vaccine
- Complete with a draw off system

PART NO	DESCRIPTION
532058	Phillips Semi. Disp. Vaccinator 5ml with Draw Off

5ML EZI-GRIP VACCINATOR WITH BOTTLE ATTACHMENT



The Ezi-Grip Vaccinator has a clear dose selection, is easy to clean and comfortable to use. Its lightweight ergonomic design with variable dosing makes it suitable for cattle, sheep and pigs. The Ezi-Grip vaccinator is constructed from high quality materials and built to last. The bottle attachment draws vaccine from the bottle or flexible pack.

KEY FEATURES INCLUDE:

- Can be sterilised without dismantling
- Luer lock needle fitting avoids accidental dislodging
- Easy to clean and comfortable to use

PART NO	DESCRIPTION
532400	5ml Ezi-Grip Vaccinator with Bottle Attachment

2ML AUTO VACCINATOR WITH BOTTLE ATTACHMENT



The Phillips 2ml Auto Vaccinator is made of metal construction, which is designed to last. This is a precision instrument, that is fully serviceable. Proven design, tested and trusted to give accurate vaccinations for years. The vaccinator comes with a bottle attachment to draw vaccine from a bottle or flexible pack. Suitable for iron injection.

KEY FEATURES INCLUDE:

- Calibration range 0.5ml to 2ml
- Spare part kits are available
- Fully serviceable
- Proven design, tested and trusted to give accurate vaccinations for years

PART NO	DESCRIPTION
530062	2ml Auto Vaccinator with Bottle Attachment

12.5ML AUTO DRENCHER/VACCINATOR



A versatile instrument, suitable for low volume drenching of sheep and injection of cattle. Very simple to use and is made from durable plastic and is fully serviceable. Provides accurate and reliable dosing.

KEY FEATURES INCLUDE:

- Simple to use
- Made from durable plastic and fully serviceable
- Provides accurate and reliable dosing
- 1ml graduations (1ml – 12.5ml)

PART NO	DESCRIPTION
440097	12.5ml Auto Drencher/Vaccinator

5ML AUTO VACCINATOR WITH BOTTLE ATTACHMENT



The Auto Vaccinator 5ml, is of a robust metal construction, which is designed to last. The Vaccinator only requires a gentle squeezing action to administer the dose. The Vaccinator comes with a bottle attachment to draw the vaccine from a bottle or flexible pack.

KEY FEATURES INCLUDE:

- Luer lock to secure needle
- Calibration 0.5ml to 5ml
- Robust metal construction, designed to last

PART NO	DESCRIPTION
531091	5ml Auto Vaccinator with Bottle Attachment

20ML AUTO DRENCHER



The 20ml Auto Drencher is one of the most popular and accurate drenchers in the world. It has a high quality, precision engineered metal construction, suitable for a variety of livestock. In-line testing during manufacture ensures the very best quality. Dosages are easily adjusted in a second for an economic and efficient use of drench.

KEY FEATURES INCLUDE:

- Quick fill for speedy operation
- Specially designed anti-spin nozzle
- Comfort grip offers convenience for the operator and ensures animals are treated with minimal fuss and effort
- Rear inlet valves for easy cleaning and suitable to be converted to an injector
- Cylinder has 1ml graduations

PART NO	DESCRIPTION
441053	20ml Auto Drencher

25ML AUTO DRENCHER



An economic, lightweight drencher made from high quality plastic for oral dosing of livestock. The Drencher is compatible with most liquid drenches.

KEY FEATURES INCLUDE:

- Compatible with most liquid drenches
- Manufactured from high grade plastic
- The 25ml Auto Drencher has 2.5ml graduations
- Easy access clip opening for cleaning
- Premium tube feed

PART NO	DESCRIPTION
439046	25ml Auto Drencher

30ML AUTO DRENCHER



Don't waste a drop of drench! Manufactured from high grade plastic for a robust construction - this drencher is virtually unbreakable. It has a bite resistant metal nozzle, adjustable for comfortable drenching. Used to administer oral animal health products to a variety of livestock. Removable cylinder makes cleaning and servicing easy.

KEY FEATURES INCLUDE:

- Removable cylinder makes cleaning and servicing easy
- Manufactured from high grade plastic for a robust construction
- The 30ml Auto Drencher has 2.5ml graduations
- Feed tube included
- Spare part kits are available

PART NO	DESCRIPTION
440048	30ml Auto Drencher

60ML PHILLIPS AUTO DRENCHER



The Phillips Auto Drencher gives an accurate dose every time, with the minimum of fuss. It is made from high grade plastic, which has rendered the drencher virtually unbreakable. It comes complete with a bite resistant metal nozzle, and is fully adjustable for comfortable drenching.

KEY FEATURES INCLUDE:

- Easy to adjust, clean and maintain
- Bite resistant metal nozzle, adjustable for comfortable drenching
- Provides accurate and reliable dosing
- 2.5ml graduations

PART NO	DESCRIPTION
440036	Phillips Auto Drencher 60ml

150ML NON AUTO CATTLE DRENCHER



Use the non Auto drencher to suck up the drench through the nozzle, set the dose and start to drench with one of the easiest and fastest drencher's around. Features a clear plastic barrel and easy dosage settings.

KEY FEATURES INCLUDE:

- Bite resistant metal nozzle
- Easy to maintain and to clean
- 25 ml graduations.

PART NO	DESCRIPTION
439034	Non Auto Cattle Drencher 150ml

MASTERJECT AUTO INJECTOR



The MasterJect cattle injection gun is an automatic, long-reach animal injection system that utilises a 20ml syringe that can be set to any dose rate from 1ml up to 20ml. The MasterJect Auto Injector allows the administration of multiple doses of medication to animals from up to 5' (1.524 metres) without having to restrain them. Its aluminium handle is easily adjusted and can be locked to the required length.

Though MasterJect is primarily a cattle injection gun, it can also medicate pigs, and any other thick-skinned animal. The beauty of MasterJect is that, unlike other injection or darting solutions, it immediately stops the flow of medication should the animal pull away from the needle.

KEY FEATURES INCLUDE:

- Long-reach injection gun
- Medicate cattle from up to 5 feet away
- Dose rates from 1ml to 20ml
- Flow stops automatically if animal pulls away from the needle
- Can also be used with pigs and other thick-skinned animals
- Each system comes complete with 4 syringes and needles and its own handy carry bag

PART NO	DESCRIPTION
AH01001	Masterject Auto Injector

RITCHEY 15ML APPLICATOR 90MM DRENCH NOZZLE




A robust drench nozzle for adult animals featuring a robust ergonomic 'V-grip' design and comes complete with draw off tube and screw thread adjuster.

KEY FEATURES INCLUDE:

- 1ml graduations from 1-15ml

PART NO	DESCRIPTION
40504	15ml Applicator 90mm Drench Nozzle

MASTERJECT NEEDLE



A pack of 12 needles used in the Masterject Automatic Animal Injector System

PART NO	DESCRIPTION
AH01003	Masterject Needles pack of 12

DISPOSABLE SYRINGES



Quality disposable syringes with clear graduations and rubber piston for excellent suction. Available individually or in packs.

PART NO	DESCRIPTION	
AH02016	Disposable Syringe 50MI	1/Each
AH02006	Disposable Syringe 50MI	Box Of 25
AH02009	Disposable Syringe 50MI	Box Of 50
AH02021	Disposable Syringe 30MI	1/Each
AH02013	Disposable Syringe 30MI	Box Of 50
AH02020	Disposable Syringe 20MI	1/Each
AH02012	Disposable Syringe 20MI	Box Of 50
AH02019	Disposable Syringe 10MI	1/Each
AH02010	Disposable Syringe 10MI	Box Of 100
AH02018	Disposable Syringe 5MI	1/Each
AH02014	Disposable Syringe 5MI	Box Of 100
AH02017	Disposable Syringe 2MI	1/Each
AH02011	Disposable Syringe 2MI	Box Of 100

6ML STV BOTTLE MOUNT INJECTOR



The Ritchey mount injector bottle fits easily and directly on to the syringe, so there is no need for tubes. It has 0.5ml graduations from 0.5ml - 6ml and a metal Luer Lock Tip to fit most Luer specification needles.

KEY FEATURES INCLUDE:

- Screw thread adjuster
- Ergonomic V-grip design

PART NO	DESCRIPTION
40502	6ml Stv Bottle Mount Injector

LUER NEEDLES PK 12



Luer Needles are used for veterinary and livestock applications. They are fully sterilisable needles and are made from high quality stainless steel.

KEY FEATURES INCLUDE:

- Fully sterilisable
- Made from high quality stainless steel

Packaged in boxes of 12

PART NO	DESCRIPTION
534241	Luer Needles 13G x 1 3/4" single
534249	Luer Needles 13G x 2" Pk 12
534234	Luer Needles 14G x 1 1/4" Pk 12
534228	Luer Needles 14G x 3/4" Pk 12
534201	Luer Needles 16G x 1 1/2" Pk 12
534195	Luer Needles 16G x 1" Pk 12
534189	Luer Needles 16G x 1/2" Pk 12
534192	Luer Needles 16G x 3/4" Pk 12
534183	Luer Needles 17G x 1" Pk 12
534158	Luer Needles 17G x 1/2" Pk 12
534171	Luer Needles 17G x 3/4" Pk 12
534146	Luer Needles 18G x 1" Pk 12
534122	Luer Needles 18G x 1/2" Pk 12
534134	Luer Needles 18G x 3/4" Pk 12
534092	Luer Needles 19G x 1/2" Pk 12
534109	Luer Needles 19G x 3/4" Pk 12
534067	Luer Needles 20G x 1/2" Pk 12
534079	Luer Needles 20G x 5/8" Pk 12
534043	Luer Needles 21G x 1/2" Pk 12
534031	Luer Needles 22G x 1/2" Pk 12
534018	Luer Needles 22G x 1/4" Pk 12

FLUTTER VALVE INJECTION SET



The Flutter Valve Injector Set allows for injection of large volumes of solution. Supplied with two record fitting needles, 13 gauge 44.6mm (1.75") - For glass or plastic bottles.

KEY FEATURES INCLUDE:

- For easy injection of large volumes of solution
- For glass or plastic bottles

PART NO	DESCRIPTION
500069	Flutter Valve Injection Set

DISPOSABLE NEEDLES PK 100



A box of 100 individually wrapped Luer lock needles.

PART NO	DESCRIPTION	
528808	Disposable Needles	16G x 1"Pk 100
528821	Disposable Needles	16G x 1 1/2"Pk 100
528766	Disposable Needles	17G x 1/2"Pk 100
528825	Disposable Needles	17G x 3/4"Pk 100
528882	Disposable Needles	18G x 1"Pk 100
528523	Disposable Needles	18G x 1 1/2"Pk 100
528754	Disposable Needles	18G x 1/2"Pk 100
528778	Disposable Needles	18G x 3/4"Pk 100
528535	Disposable Needles	19G x 1"Pk 100
528791	Disposable Needles	20G x 1/2"Pk 100

60ML AUTO DRENCHER C/W FLOATING HOOK NOZZLE



The Phillips 60ml Auto Cattle Doser with feed tube provides an accurate dose every time with the minimum of fuss. Made using high grade plastic construction, this drencher is extremely tough and features a bite resistant metal nozzle which is fully adjustable for comfortable drenching.

PART NO	DESCRIPTION
40509	60ml Auto Drencher C/W Floating Hook Nozzle

RITCHEY FOOTCARE & CLEANSING SPRAY*

A general purpose, alcohol based, purple cleansing spray for use on all animals. Navels should be sprayed as soon as possible after birth and allowed to dry before turnout.

KEY FEATURES INCLUDE:

- For use on all animals
- Alcohol based, purple cleansing spray
- 500ml container



PART NO	DESCRIPTION
41141	Ritchey Footcare and Cleansing Spray 500 ML

* Sold in packs of 12

RITCHEY HOOFTREAT PLUS

Made using a liquid formula which helps to relieve the symptoms of many cattle lameness diseases, particularly Digital Dermatitis. Ritchey Hoof Treat Plus is tough, durable and easily used to clean hooves and surrounding areas. It also contains a dye to aid visibility on the hoof, meaning it is the new generation of non-antibiotic biocide which is free from copper sulphate and other harmful organic acids.

KEY FEATURES INCLUDE:

- Non-antibiotic biocide
- Targets & prevents many illnesses
- Aids visibility by dyeing the areas it is applied too
- Free from copper sulphate and other harmful organic acids

PART NO	DESCRIPTION
M41004	Ritchey Hoof Treat Plus 25Lt - Not illustrated

TRAFALGAR LIGHTWEIGHT FOOTROT SHEARS

The Trafalgar Footrot Shears are high quality and lightweight. They have a strong cut, are easy to handle and manoeuvre around cleats.

KEY FEATURES INCLUDE:

- Includes wrist strap
- High quality light weight shears

PART NO	DESCRIPTION
421352	Trafalgar Lightweight Footrot Shears




HOOF CUTTER

The Hoof Cutter comes complete with fitting screw and exchangeable blades and is built with double-action compound joints.

KEY FEATURES INCLUDE:

- Soft grip handles
- Designed for heavy-duty use
- Interchangeable cutting blades
- Length: 41cm

PART NO	DESCRIPTION
80010	Hoof Cutter



TRAFALGAR HEAVY DUTY FOOTROT SHEARS (SERRATED)

The Trafalgar Heavy Duty Footrot Shears are made from high quality tempered steel. They have heat treated precision ground blades for strength and supreme sharpness.

KEY FEATURES INCLUDE:

- Serrated blades
- High quality tempered steel

PART NO	DESCRIPTION
421261	Trafalgar Heavy Duty Footrot Shears (Serrated)



AESCULAP HOOF KNIVES

The Aesculap Hoof Knives are German engineered knives with an extremely effective cutting edge and are famous for their resilience and resistance to damage. The ergonomic handle allows the knives to be used comfortably even after prolonged periods of time.

Manufactured from special break-proof stainless steel which gives the blade its long life. Making these one of the most popular ranges of cattle and horses hoof knives available today.

KEY FEATURES INCLUDE:

- Ergonomic handle
- Long life stainless steel blade
- Hard-wearing 'break-proof' Aesculap hoof knife
- Appreciated by professional users for many years
- Length 19cm
- Left and right handed knives available

PART NO	DESCRIPTION
509012	Hoof Knife R/H Single Edge
509024	Hoof Knife L/H Single Edge
509048	Hoof Knife Swedish Style R/H Single Edge
509061	Hoof Knife Swedish Style L/H Single Edge
509073	Hoof Knife Swedish Style R/H Double Edge
509115	Hoof Knife Small Double Edge
509139	Hoof Knife R/H Double Edge
509176	Hoof Knife Aesculap L/H Single Edge Narrow
509188	Hoof Knife Aesculap R/H Double Edge
509206	Hoof Knife Aesculap R/H Single Edge Narrow
509210	Hoof Knife Aesculap R/H Double Edge Medium
509212	Hoof Knife Aesculap R/H Single Edge Narrow
509218	Hoof Knife Aesculap L/H Single Edge Narrow



PROFI HOOF KNIVES

The Profi Hoof Knife is made of high quality special steel developed for cutting tools. The blade is easy to sharpen, and includes an improved blade radius for precise action. The ergonomic shaped handles are made of fine wood.

KEY FEATURES INCLUDE:

- High quality steel, for a longer lasting edge and easier sharpening
- Broad blade and fine wooden handle
- Strong and hard wearing

PART NO	DESCRIPTION
80001	PROFI Hoof Knife Right Handed Thin Blade
80002	PROFI Hoof Knife Left Handed Thin Blade
80003	PROFI Hoof Knife Double Edged Medium Blade



Lambing



COX EWE TRUSS*

The Cox Ewe Truss prevents protruding or forced bodies during the lambing process. It is made from rot resistant nylon webbing, with Kevlar stitching and with quick release buckles. It allows ewes to lamb normally.

KEY FEATURES INCLUDE:

- Prevents protruding or forced bodies
- Made from rot resistant nylon webbing, with Kevlar stitching
- Quick release buckles make it safe and easy to use



PART NO	DESCRIPTION
434073	Cox Ewe Truss

* Sold as a single item

RITCHEY EWE TRUSS*

The Ritchey Ewe Truss can help to prevent a uterine prolapse and will allow lambing to proceed, whilst exerting pressure on key points. It can also be used to hold back and support a prolapsed uterus, which will lessen the risk of infection and speed up the recovery process. The device keeps the head raised and back down to relieve pain.

KEY FEATURES INCLUDE:

- Made from washable, man-made material that won't rot
- Eases the lambing process
- Exerts pressure on key points
- Keeps head raised and back down to relieve pain
- Prevents the ewe's uterus from prolapsing



PART NO	DESCRIPTION
58000	Ritchey Ewe Truss

* Sold as a single item

UTERI SUPPORT*



The Uteri Support is an internal uterine support system to help prevent ewe prolapses during birth. The product can be used on its own, by tying it to the fleece, or alternatively in conjunction with an ewe truss. Always sterilise the product before use, and apply Ritchey Superlube before insertion.

KEY FEATURES INCLUDE:

- Eases pain when giving birth
- Reduces the risk of a uterine prolapse

PART NO	DESCRIPTION
BR04070	Uteri Support

* Uteri support sold as a pack of 4

LAMBING AID SET



The Lambing Aid set is a useful aid during a difficult lambing. The set includes a 3mm (1/8") lambing rope and sensory fingertips.

KEY FEATURES INCLUDE:

- 3mm (1/8") lambing rope is included
- Sensory fingertips included

PART NO	DESCRIPTION
434085	Lambing Aid Set

* Sold as a single item

LAMBING ROPE*



The Lambing Rope is designed to provide safe, gentle assistance during birth. The Lambing Rope is made from nylon and has a single loop on the end. The product is also durable enough to withstand cleaning with disinfectant products, essential for lambing season.

KEY FEATURES INCLUDE:

- Nylon with a single loop
- Diameter 3mm (1/8")
- Can be washed and reused

PART NO	DESCRIPTION
434024	Lambing Rope

* Sold as a single item

ELASTRATOR LAMB CASTRATION RINGS



The Elastrator Lamb Castration rings are made from natural latex, and allow you to dock the lamb's tail. Docking improves the health and welfare of sheep and lambs, preventing faecal matter from accumulating on the tail and hindquarters of the animal. Apply a rubber ring to the tail using the elastrator tool, the bands will then cut off the blood supply to the tail, causing the tail to fall off in 7 to 10 days.

KEY FEATURES INCLUDE:

- Rings are made from natural latex
- Quality you can rely on
- Trusted the world over for over 50 years
- Can prevent faecal matter accumulating on the tail

PART NO	DESCRIPTION
426027	Elastrator Rings Box 100
426039	Elastrator Rings Box 500
426118	Elastrator Rings Applicator

LAMBING INSTRUMENT*



The Lambing Instrument is vital to have around the farm when an ewe is experiencing a difficult birth, and especially when the lamb's head is back, the lambing instrument can help to provide some valuable extra reach. The loop at the top is designed to fit over the lamb's head, whilst the fork-shape supports the animal's chin.

KEY FEATURES INCLUDE:

- Lightweight, durable and flexible
- Extremely easy to use
- Purpose-built to help with lambing, and an essential aid

PART NO	DESCRIPTION
434036	Lambing Instrument

* Sold as a single item

CASTRATION RINGS*



Castration and tail docking should be carried out within the first week of life, and care must be taken to place the ring above both testes.

KEY FEATURES INCLUDE:

- Made from rolled rubber
- Available in two different quantities

PART NO	DESCRIPTION
56402	Castration Rings 500
56404	Castration Rings Bucket 1500

* Castration rings are sold in packs of 10

CASTRATION RING PLIERS



The Castration Ring Pliers are made with the highest quality materials and are strong and reliable with good length of prongs. A Patented recess helps to avoid the rings from slipping off.

KEY FEATURES INCLUDE:

- High quality pliers
- Strong and reliable patented recess to keep rings secure
- Complete with 100 rings

PART NO	DESCRIPTION
56411	Castration Ring Pliers

FARMER'S CHOICE LAMB COLOSTRUM*

FARMER'S CHOICE



Specifically formulated to provide lambs with the best possible start. The Colostrum provides four instant energy sources in the form of carbohydrates. Carbohydrates, unlike most fats, are easily converted into glucose and energy to boost lambs and get them on their feet to suckle.

KEY FEATURES INCLUDE:

- Farmer's Choice has been developed by farmers, for farmers and has been trialled on farms across the UK
- Colostrum contains natural nucleotides and, using breakthrough technology, it draws on advances in human baby nutrition

PART NO	DESCRIPTION
710010	Lamb Colostrum (10 feed) - 250 gsm
710020	Lamb Colostrum (20 Feed) - 500 gsm
710050	Lamb Colostrum (100 Feed) 2.5 kg
710060	Lamb Colostrum (400 Feed) 10 kg

* 500 gsm - Sold in packs of 6

EMERGENCY FEEDER



The Emergency Feeder is designed to be used for emergencies only, and provides a new born lamb or calf with the colostrum or glucose they need, when they are unable to get it any other way. This device is designed to go down into a newborn's stomach, and there are safeguards in place to prevent it from entering the trachea (windpipe).

KEY FEATURES INCLUDE:

- Suitable for lambs or calves
- Easy to fill, clean and use
- Has an extra-large capacity
- Leads directly to the stomach
- 140ml Syringe
- Supplied with two catheters

PART NO	DESCRIPTION
56011	Emergency Feeder – 140ml Lamb

* Sold as a single item



MADE IN GREAT BRITAIN

71000
71001

71132
71130

71230

71032
71030

71022
71020

ROB-ROY CROOKS

Available with a long or short neck, the shaft of the Rob Roy Crook is made from an aluminium poly coated tube.

KEY FEATURES INCLUDE:

- Available in 2 lengths
- Strong, durable construction
- Comfortable ergonomic handle

Available in red

PART NO	DESCRIPTION
71020	Crook: Rob-Roy Long Neck
71022	Crook: Rob-Roy Short Neck

Sold as a single item

BLACK COL-ROY CROOKS

Col-Roy Crooks are made with a unique fibre-glass flexible shaft and a strong nylon head and are available in two sizes depending on requirements.

KEY FEATURES INCLUDE:

- Available in 2 lengths (long 134cm and medium 118cm)
- Made from unique fibre glass
- Flexible Shaft
- Strong nylon head

PART NO	DESCRIPTION
71000	Crook: Col-Roy Original Long
71001	Crook: Col-Roy Medium

Sold as a single item

BLU-ROY CROOKS

The Blu-Roy Crook is made from solid aluminium with cobalt blue shaft and is polycoated for comfort. It is available in 5 different sizes depending on requirements.


KEY FEATURES INCLUDE:

- Available in 5 sizes
- Strong, durable construction
- Comfortable ergonomic handle

PART NO	DESCRIPTION
71030	Crook: Blu-Roy Long Neck
71032	Crook: Blu-Roy Short Neck
71130	Crook: Blu-Roy Long Leg
71132	Crook: Blu-Roy Short Leg
71230	Crook: Blu-Roy Combination

Sold as a single item

SYRINGE WITH CATHETER



60ml Syringe and Catheter, which is used for feeding lambs that are too weak to suckle due to exposure or other ailments.

KEY FEATURES INCLUDE:

- 60ml syringe
- Helps to feed lambs which are otherwise too weak

PART NO	DESCRIPTION
435107	Syringe with catheter (Lamb)
435070	Syringe with rubber catheter

STERIMATIC REFILL



A refill pack for use with the Disposable Vaccinator.

KEY FEATURES INCLUDE:

- Five Stericaps and needles are included

PART NO	DESCRIPTION
530190	Sterimatic Refill

COX 5ML STERIVAX PACK



The Sterivax is a precision moulded plastic container that fits into the Sterimatic sleeve and contains foam soaked in disinfectant, which cleans the needle as it passes through the Stericap.

KEY FEATURES INCLUDE:

- The Sterivax pack contains: 5ml semi-disposable vaccinator, draw off system, 1 x Sterimatic Sleeve, 3 x Stericaps and 3 x 18g x 1" needles.
- Give up to 300 automatically disinfected and accurate injections
- Increased safety for the operator allowing for fewer needle stick and scratch injuries
- Healthier livestock because of reduced abscesses and cross infections between animals

PART NO	DESCRIPTION
530311	Cox 5ml Sterivax Pack

SEMI DISPOSABLE VACCINATOR 5ML



A High quality precision Phillips instrument. The Semi Disposable Vaccinator is one of the most popular automatic plastic vaccinators in the UK. It is lightweight and economic and provides accurate administration of vaccine. This injector is used as part of the popular Sterivax system.

KEY FEATURES INCLUDE:

- High quality precision vaccinator
- Lightweight and economic
- Provides accurate administration of vaccine
- Comes complete with draw off system

PART NO	DESCRIPTION
532058	Semi Disposable Vaccinator 5ml

FOR MORE
DRENCHERS
& VACCINATORS
GO TO
PAGE 15

450ML LAMB BOTTLE

Always needed around the farm. Supplied with teat and featuring a large bottle neck for easy filling, as well as with clear graduations for accurate dosage.



KEY FEATURES INCLUDE:

- Complete with teat
- Individually adjustable milk flow
- Bottle neck with large diameter for easy filling
- Large graduations
- 450ml capacity

PART NO	DESCRIPTION
435032	450ml Lamb Bottle

* Sold as a single item

LAMB FEEDER BOTTLE WITH TEAT

The Lamb Feeder Bottle is a must have piece of equipment come lambing season, the Ritchey Lamb Feeding Bottle comes complete with a soft, screw-on, non-vac teat, and is the perfect tool for hand rearing new born lambs.



KEY FEATURES INCLUDE:

- Made from a durable plastic
- Red teat can be changed and replaced if it gets too worn or damaged
- Includes an in-built flutter valve to prevent any liquid or air making its way back in to the bottle
- Incorporates a natural feeling soft, red teat

PART NO	DESCRIPTION
56230	Lamb Feeder Bottle with Teat
56207	Spare Teat – 5 per pack

* Sold as a single item

LAMB BUCKET FEEDER AND 6 TEATS

Lamb Feeder System and hanging rack, incorporating a traditional cross cut teat system to eliminate wastage. It is a square plastic feeder with 2 chambers, which hangs easily on a rail.



KEY FEATURES INCLUDE:

- 10 litre capacity
- Includes carry handle and 6 teats
- Can be fixed onto fences or feeding rails

PART NO	DESCRIPTION
503241	Lamb Bucket Feeder and 6 Teats
503253	Spare teat

* Sold as a single item

OVER HURDLE 6 TEAT LAMB FEEDER

Square bucket, complete with 6 teats with non-return valves to ensure a steady milk supply. Teat holes can be enlarged for greater milk flow. For positioning on a pen partition.



KEY FEATURES INCLUDE:

- Plastic carry handle
- For feeding up to 6 lambs
- Holds up to 12 litres

PART NO	DESCRIPTION
503245	Special 6 Teat Lamb Feeder
56128	Spare Valve Kits
56133	Valve Washers

* Sold as a single item

3 TEAT LAMB FEEDING BUCKET

Complete with 3 teats.

PART NO	DESCRIPTION
503200	3 Teat Lamb Feeding Bucket
505134	Spare Teats



LAMB BAR FEEDER WITH NON-RETURN VALVE



The Lamb Bar Feeder is a functional feeding system, that can be attached to a fence or feeding rail. The product comes complete with a teat, tubing, and a non-return valve to ensure a steady milk supply.

KEY FEATURES INCLUDE:

- Fitted with non-return valve
- Can be fixed onto fences or feeding rails

PART NO	DESCRIPTION
504142	Lamb Bar Feeder with N R Valve

BLACK PULL ON LAMB TEAT



Made of superior rubber, making it supple yet strong. A bulbous end, which is specially designed to encourage fast acceptance. Hollow cone tip permits milk flow regulation, by cutting the tip for a slow flow or at the lower tip for increased flow.

KEY FEATURES INCLUDE:

- Fits most household bottles
- Bulbous end to encourage fast acceptance

PART NO	DESCRIPTION
503101	Black Pull On Lamb Teat

OVER HURDLE SPARE TEATS

PART NO	DESCRIPTION
56203	Spare Teats Red - Not Illustrated
56204	Spare Teats White - Not Illustrated

LAMB TEAT PLUG IN



The Lamb Teat Plug In is pushed inside the neck of a bottle to assist with lamb feeding. The milk hole can be enlarged to your requirements.

KEY FEATURES INCLUDE:

- Rubber teat, which is pushed inside the neck of a bottle
- Milk hole can be enlarged where necessary

PART NO	DESCRIPTION
503137	Lamb Teat Plug In

PRITCHARDS SCREW ON LAMB TEAT

A flutter valve lets air back into the bottle as the lamb is drinking. Seals onto most modern bottles.

KEY FEATURES INCLUDE:

- Rubber teats for bottles with a threaded neck
- Fits many plastic 'pop' bottles



PART NO	DESCRIPTION
503228	Pritchards Screen On Lamb Teat

SCREW ON LAMB TEAT *

Soft silicone rubber and metal ball valve incorporated in the cap.

KEY FEATURES INCLUDE:

- Fits most plastic bottles



PART NO	DESCRIPTION
56207	Screw on Lamb Teat

* Teats sold in packs of 5

COX SCREW ON TEAT



The Cox Screw on teat fits the Feeding Bottle and most other plastic bottles and is ideal for Alpacas, Goats & Lambs.

KEY FEATURES INCLUDE:

- Latex teats for bottles with a threaded neck
- Fits most plastic 'pop' bottles

PART NO	DESCRIPTION
503232	Cox Screw On Teat

INFRARED REFLECTOR



With 8 ventilation holes to provide better air circulation, the aluminium shade helps to extend the lamps service life; whilst the attachment chain allows it to hang securely.

KEY FEATURES INCLUDE:

- Aluminium shade, diameter 21cm
- Eight ventilation holes provide better air circulation, extending the lamp's service life
- Chain for attachment, max load 20kg
- Comes with 5m cable

PART NO	DESCRIPTION
56810	Infrared Reflector

PHILIPS INFRARED ECONOMY LAMP



PART NO	DESCRIPTION
56813	Philips Infrared Economy Lamp Red 100W
56814	Philips Infrared Economy Lamp Red 175W

POLYTHENE WRIST LENGTH GLOVES



Polythene, wrist length gloves (30cm). One size fits all. Comes in a dispenser bag of 100 gloves. Always useful to have around.

KEY FEATURES INCLUDE:

- Inexpensive clear polythene gloves, for general hygiene on the farm
- Available in packs of 100

FOR MORE GLOVES GO TO PAGE 5

PART NO	DESCRIPTION
496157	Polythene Wrist Length Gloves

INFRARED LAMP RED



To be used with the Infrared Reflector.

KEY FEATURES INCLUDE:

- Hard glass
- Tolerates splash water

Available in 150W and 250W

PART NO	DESCRIPTION
56811	Infrared Lamp Red 150W
56812	Infrared Lamp Red 250W

LAMB COVERS



The polythene lamb jackets help prevent hypothermia, which is a major cause of lamb mortality in the UK. The Lamb jackets are designed not to interfere with mothering, as the ewe will still be able to smell and identify her own lamb, due to the cut away back.

KEY FEATURES INCLUDE:

- Polythene lamb jacket which protects new born lambs against the elements
- Includes a marking pen for numbering lambs
- Cut away back ensures no interference with the natural bonding process
- As the lamb grows the need for protection decreases and it can burst through and shed the biodegradable cover
- 100 lamb jackets per box
- Environmentally friendly - the polythene jacket is degradable in ultra-violet light

PART NO	DESCRIPTION
56802	Lamb Covers (100) Single Orange
56803	Lamb Covers (100) Twin Clear

PLEASE NOTE: When purchasing our lamb jackets please bear in mind that lambs will typically shed their jacket after 3-days, so please adjust your order quantity accordingly.

GAMBREL RESTRAINER



The Gambrel Restrainer is the easy and simple way to restrain and hold sheep. It is an essential aid for lambing, holding, foot trimming, treatment and transportation. For use with growing lambs and large breeds of sheep.

KEY FEATURES INCLUDE:

- Easy to use and quick to fit
- 10% holding, 90% tranquillity, it keeps them so tranquil they will even feed
- Easy to use, humane way to restrain stock

PART NO	DESCRIPTION
437177	Gambrel Restrainer Large

RITCHEY SUPER SPRAYLINE



Ritchey Super Sprayline is an all-weather, long lasting aerosol marker for general use in sheep and lamb identification. It is also used as a multi-purpose dairy cow flank and udder marker. Super Sprayline is approved by the British Wool Marketing Board.

KEY FEATURES INCLUDE:

- Can be used on cows, goats, sheep and pigs
- Long term spray marker
- 360 degree spray
- Suitable for use on wet or dry fleeces
- Block free nozzle
- Ideal for sorting, grading, scanning or pairing ewes with lambs after birth
- Available in a 200ml or 400ml can
- Colours Purple, Yellow, Red, Orange, Blue, Green and White (White available in 400ml only)



PART NO	DESCRIPTION
30003	Super Sprayline 200ml Purple
30004	Super Sprayline 200ml Yellow
30005	Super Sprayline 200ml Red
30006	Super Sprayline 200ml Orange
30007	Super Sprayline 200ml Blue
30008	Super Sprayline 200ml Green
30102	Super Sprayline 400ml White
30103	Super Sprayline 400ml Purple
30104	Super Sprayline 400ml Yellow
30105	Super Sprayline 400ml Red
30106	Super Sprayline 400ml Orange
30107	Super Sprayline 400ml Blue
30108	Super Sprayline 400ml Green

PLEASE NOTE: 200 ml Sprayline is sold in packs of 20
400 ml Sprayline is sold in packs of 12

RITCHEY STOCK MARKER SPRAY 450ML



Available in 3 vividly clear colours, Ritchey Stock Marker Spray can be used on cows and sheep, etc. and is designed to be sprayed directly onto an animal so that it can be clearly identified and singled out when necessary.

KEY FEATURES INCLUDE:

- Great value for money
- Can be used on cows or sheep
- Creates durable short-term markings
- Available in Red, Blue and Green

PART NO	DESCRIPTION	
30315	Ritchey Stock Marker Spray 450ml Red	
30317	Ritchey Stock Marker Spray 450ml Blue	
30318	Ritchey Stock Marker Spray 450ml Green	

DALES MARKER



The Dales Sheep Marker has been designed to cater for upland sheep, which are generally more prone to experiencing all elements of the weather. Dales Sheep marker fluid adheres strongly to a sheep's coat and is available in a range of colours, making it ideal for numbering/ marking your flock.

KEY FEATURES INCLUDE:

- Versatile and durable fleece marker
- Able to adhere directly onto sheep whether wet or dry
- Can be used pre or post dipping
- Can be used alongside branding irons & spot marking sticks
- Available in 2.5lt and a 5lt tub
- Available in Black, Purple, Orange, Blue, Green and Deep Red.

PART NO	DESCRIPTION	
33101	Dales Marker 2.5 Lt Black - sold in packs of 4	
33103	Dales Marker 2.5 Lt Purple - sold in packs of 4	
33106	Dales Marker 2.5 Lt Orange - sold in packs of 4	
33107	Dales Marker 2.5 Lt Blue - sold in packs of 4	
33108	Dales Marker 2.5 Lt Green - sold in packs of 4	
33112	Dales Marker 2.5 Lt Deep Red - sold in packs of 4	
33201	Dales Marker 5 Lt Black - sold in packs of 2	
33203	Dales Marker 5 Lt Purple - sold in packs of 2	
33204	Dales Marker 5 Lt Yellow - sold in packs of 2	
33206	Dales Marker 5 Lt Orange - sold in packs of 2	
33207	Dales Marker 5 Lt Blue - sold in packs of 2	
33208	Dales Marker 5 Lt Green - sold in packs of 2	
33212	Dales Marker 5 Lt Deep Red - sold in packs of 2	

PLEASE NOTE:

- 2.5lt sold in cases of 4
- 5lt sold in packs of 2



COX SHEEP MARKER SPRAY



Cox Sheep Marker Spray is Ideal for grading or pairing ewes and lambs at birth.

KEY FEATURES INCLUDE:

- Can be applied wet or dry to wool
- Will remain bright even after heavy rain
- Specifically designed tip for precise and accurate markings
- High quality – accredited and approved by the British Wool Marketing Board
- Available in 400ml cans
- Available in Red, Blue, Green, Purple, Black, Yellow and Orange

PART NO	DESCRIPTION	
476109RED	Cox Sheep Market Spray RED 400ml	
476109BLU	Cox Sheep Market Spray BLUE 400ml	
476109GRE	Cox Sheep Market Spray GREEN 400ml	
476109PUR	Cox Sheep Market Spray PURPLE 400ml	
476109BLA	Cox Sheep Market Spray BLACK 400ml	
476109YEL	Cox Sheep Market Spray YELLOW 400ml	
476109ORA	Cox Sheep Market Spray ORANGE 400ml	

PLEASE NOTE: Black is not BWMB approved.

WAX MARKER STICK



Wax marker sticks with twist-up action making it ideal for temporary identification.

KEY FEATURES INCLUDE:

- Resealable cap preventing accidental marks and stops the stick from becoming dried out
- A reliable and economical method of identification
- A handy pocket-sized marker
- Available in Black, Blue, Green, Purple and Red

PART NO	DESCRIPTION	
475609RED	WAX MARKER STICK RED	
475609BLU	WAX MARKER STICK BLUE	
475609GRE	WAX MARKER STICK GREEN	
475609PUR	WAX MARKER STICK PURPLE	

NYLON SUTURE TAPE




A nylon suture tape on card for use in closing a tissue opening used with suture needles. Often used for stitching a prolapse in sheep.

KEY FEATURES INCLUDE:

- Length 5m and available in 3mm or 6mm width

PART NO	DESCRIPTION
491317	Nylon Suture Tape 3mm
491305	Nylon Suture Tape 6mm

SUTURE NEEDLE HALF CURVE 6" PK 2



A high grade stainless steel suture needle which ensures a smooth, less traumatic passage through tissue.

KEY FEATURES INCLUDE:

- Stainless Steel
- 154mm (6") half curved

Available in packs of two

PART NO	DESCRIPTION
491329	Half Curve Suture Needle 6"

DIGITAL THERMOMETER




Provides a rapid and accurate reading of the temperature of all animals. Especially useful for diagnosing hypothermia in lambs.

KEY FEATURES INCLUDE:

- Digital dial making it easier to read
- Quick and accurate results
- Available in Celsius (digital) and Fahrenheit (digital)

PART NO	DESCRIPTION
520238	Celsius (digital)
520251	Fahrenheit (digital)

MEDICLAMP



The Mediclamp is used to accelerate the drying process of the umbilical cord. This then falls off within a shorter time frame and grants the new-born animal perfect protection to avoid bacterial infection.

As soon as the umbilical cord has fallen off, the Mediclamp seals the other end with permanent pressure.

KEY FEATURES INCLUDE:

- The Mediclamp is made of a high quality special plastic material and is designed for a long life and high reusability
- Snap on at birth to help avoid deadly diseases such as e.coli and joint ill
- Reduces the need for sprays, dips and expensive antibiotics
- Fully re-usable after sterilisation
- Available in packs of 50

PART NO	DESCRIPTION
54201	Mediclamp

ONLY THE BRIGHTEST EWES GET TOP MARKS

USE TOP MARKS HARNESS AND CRAYONS DURING MATING FOR CLEARER, BRIGHTER AND LONGER LASTING RESULTS



The Top Marks Harness has been specifically designed for comfort and ease of use. The harness moulds to the body of the ram and thus helps prevent chafing whilst securely holding the crayon firmly in place. The Top Marks Harness helps to protect your investment in your breeding stock whilst providing an accurate mark.



HARNESS KEY FEATURES:



- ✓ Comfortable Fit
- ✓ Hugs the ram's brisket and holds the crayon in place therefore reducing chafing
- ✓ Uses Quick Clips instead of buckles, making it faster and easier to strap on
- ✓ Harnesses can be strapped on without tipping the ram
- ✓ Made with rot proof Kevlar stitching
- ✓ Bright webbing makes rams easier to spot in the field

CRAYON KEY FEATURES:



- ✓ More effective than raddle paint
- ✓ Leaves a clear mark
- ✓ Helps to provide a greater breeding insight
- ✓ Click-in-Click-out Crayons do not require a cotter pin
- ✓ Bright colours and easy access packaging

TOP MARKS HARNESS & RAM CRAYONS*		
PART NO	DESCRIPTION	
477072	Top Marks Harness	
477057BLU	Top Marks Crayon All Weather Blue	
477057GRE	Top Marks Crayon All Weather Green	
477057ORA	Top Marks Crayon All Weather Orange	
477057PUR	Top Marks Crayon All Weather Purple	
477057RED	Top Marks Crayon All Weather Red	
477057YEL	Top Marks Crayon All Weather Yellow	

TOP MARKS COLD WEATHER CRAYONS		
PART NO	DESCRIPTION	
477059BLU	Top Marks Crayon Cold Blue	
477059GRE	Top Marks Crayon Cold Green	
477059ORA	Top Marks Crayon Cold Orange	
477059PUR	Top Marks Crayon Cold Purple	
477059RED	Top Marks Crayon Cold Red	
477059YEL	Top Marks Crayon Cold Yellow	

* Top Marks Ram Crayons are sold in packs of 10 per colour

DONAGHY'S HARNESS & RAM CRAYONS

Donaghy's Ram Crayons are a recognised market leader. Manufactured in New Zealand, they have provided generations of farmers with outstanding performing Ram Crayon used in the field, a tradition which continues to this day.



KEY FEATURES INCLUDE:

- Clip-in-system for rapid crayon change
- Non-Slip adjustable buckles and V neck fitting for comfort
- Ideal for marking ewes as they are served
- All Donaghy's crayons are approved by the British Wool Marketing Board

DONAGHY'S HARNESS

PART NO	DESCRIPTION
50200	Donaghy's Harness

RAM CRAYONS SOFT* ❄️

PART NO	DESCRIPTION	
50005	Ram Crayon SOFT - Donaghys Red	🔴
50006	Ram Crayon SOFT - Donaghys Orange	🟠
50007	Ram Crayon SOFT - Donaghys Blue	🔵
50008	Ram Crayon SOFT - Donaghys Green	🟢

RAM CRAYONS HARD* ☀️

PART NO	DESCRIPTION	
50105	Ram Crayon HARD - Donaghys Red	🔴
50106	Ram Crayon HARD - Donaghys Orange	🟠
50107	Ram Crayon HARD - Donaghys Blue	🔵
50108	Ram Crayon HARD - Donaghys Green	🟢

RAM CRAYONS KEY:

- ❄️ To be used when the temperature is below 16°
- ☀️ To be used when the temperature is above 16°

* Donaghys Ram Crayons are sold in packs of 5 per colour

PLEASE NOTE: Crayons should be changed every three weeks to easily identify served ewes.

RITCHEY RADDLE POWDER



A long lasting coloured paste which is ready to mix and is easy to apply to the ram's chest for marking ewes at mating time. Raddle Powder is used as an alternative to Ram Crayons.

KEY FEATURES INCLUDE:

- Will not detach on the ram's chest or between applications
- Easy to apply
- 1 Lt of vegetable oil will blend 2.2kg of raddle powder
- Available in four colours: Yellow, Red, Blue and Green

PART NO DESCRIPTION

50524	Ritchey Raddle Powder 450g Yellow	🟡
50525	Ritchey Raddle Powder 450g Red	🔴
50527	Ritchey Raddle Powder 450g Blue	🔵
50528	Ritchey Raddle Powder 450g Green	🟢
50504	Ritchey Raddle Powder 1Kg Yellow	🟡
50505	Ritchey Raddle Powder 1Kg Red	🔴
50507	Ritchey Raddle Powder 1Kg Blue	🔵
50508	Ritchey Raddle Powder 1Kg Green	🟢
50515	Ritchey Raddle Powder 1.5Kg Red	🔴
50517	Ritchey Raddle Powder 1.5Kg Blue	🔵
50518	Ritchey Raddle Powder 1.5Kg Green	🟢

ALERT TAIL PAINT*

Alert tail paint is a unique formulation to improve mating management including heat detection without irritating the cows skin.



KEY FEATURES INCLUDE:

- Paint will be broken or rubbed off after the cows are ridden, indicating that they have cycled
- Outstanding performance in both wet and dry conditions lasting up to 20 days
- Easy to apply
- Depending on usage, coverage is between 30-40 cows with a 15cm line
- Available in five colours: Yellow, Red, Orange, Blue & Green

PART NO DESCRIPTION

50304	Alert Tail Paint 500ml Yellow	🟡
50305	Alert Tail Paint 500ml Red	🔴
50306	Alert Tail Paint 500ml Orange	🟠
50307	Alert Tail Paint 500ml Blue	🔵
50308	Alert Tail Paint 500ml Green	🟢

* Alert Tail Paint is sold in packs of 12 per colour

KAMAR HEATMOUNT DETECTORS

Kamar Heatmount Detectors are a valuable heat detection aid which can assist to identify cows ready for artificial breeding. Standing heat behaviour is the most reliable sign that a cow is ready to be bred. The detectors are a pressure sensitive device with a built-in timing mechanism designed to be activated by standing heat behaviour.



Glued onto the sacrum (tail head), pressure from the brisket of a mounting animal requires approximately 3 seconds to turn the detector from white to red. This timing mechanism helps distinguish between true standing heat versus false mounting activity.

KEY FEATURES INCLUDE:

- Detectors provide a clear indication of a cow in heat by turning bright red, clearly visible from a distance
- Can reduce the cost of missed heats
- Stays firmly glued to the cow
- Helps to identify night heats
- Reliable even in extreme weather
- Increased herd profitability

PART NO	DESCRIPTION	PACK QTY	
485081	Kamar Red Heatmount Detectors	25	
485093	Kamar Blue Heatmount Detectors	25	

TAIL TAPE 50M*



KEY FEATURES INCLUDE:

- Instant visual aid for cows under observation
- Roll length 50m (165ft)
- Available in six colours: Black, Blue, Green, Red, White and Yellow

PART NO	DESCRIPTION	LENGTH	
473029BLA	Tail Tape Plain Black	50m (165ft)	
473029BLU	Tail Tape Plain Blue	50m (165ft)	
473029GRE	Tail Tape Plain Green	50m (165ft)	
473029RED	Tail Tape Plain Red	50m (165ft)	
473029WHI	Tail Tape Plain White	50m (165ft)	
473029YEL	Tail Tape Plain Yellow	50m (165ft)	

* Sold as a single item

TAIL PAINT BRUSH ON*

Donaghy's uniquely formulated liquid tail paint has been used for years as a reliable tool for heat detection in cows and a practical solution for managing herd mating.



KEY FEATURES INCLUDE:

- Oil based formula lasts up to 30 days in suitable conditions
- Trials have shown it lasts well and is easy to apply even in wet weather
- Available in five bright non-fading colours; Yellow, Red, Blue, Green and Pink

PART NO	DESCRIPTION	
50314	Tail Paint Brush On 1L Yellow	
50315	Tail Paint Brush On 1L Red	
50317	Tail Paint Brush On 1L Blue	
50318	Tail Paint Brush On 1L Green	
50319	Tail Paint Brush On 1L Pink	

- * The 1Lt Bottle is supplied with handy brush applicator cap
- * Tail Paint Brush On is sold in packs of 12 per colour

SI-RO-MARK

THE WORLDS PREMIUM SHEEP BRANDING FLUID



DURABLE	Stays on the fleece for up to one year
SCOURABLE	Approved by the British Wool Marketing Board
ECONOMICAL	Easy to use and can brand up to 1,000 sheep from one 5 litre can
DISTINCT	Available in red, blue, orange, green, purple and black
FLOCK SECURITY	Clear deterrent against theft



PART NO	DESCRIPTION	
475440BLU	Si-Ro Mark Sheep Fluid BLUE box of 4 x 5lt	
475440GRE	Si-Ro-Mark Sheep Fluid GREEN Box of 4 x 5lt	
475440ORA	Si-Ro-Mark Sheep Fluid ORANGE Box of 4 x 5lt	
475440PUR	Si-Ro-Mark Sheep Fluid PURPLE Box of 4 x 5lt	
475440RED	Si-Ro-Mark Sheep Fluid RED Box of 4 x 5lt	
475446BLA	Sheep Brand Fluid BLACK 4 x 5lt	

SHEPHERD'S WHISTLES



PART NO	DESCRIPTION
552112	Mouth Whistle Plastic
552120	Mouth Whistle Metal - Not illustrated
552094	Mouth Whistle Standard Pitch Dog - Not illustrated
552100	Mouth Silent Dog Whistle - Not illustrated

S.W.A.T AUTO FLY CONTROL SYSTEM



S.W.A.T. is an automatic fly control system, consisting of an automatic dispenser and aerosol Fly Killer. Battery powered, the unit can be set up anywhere on Farm where needed, it will work for over a month on a single refill, providing complete peace of mind for the user, with a powerful, yet natural Pyrethrin based broad spectrum insecticide which is, non-toxic and safe to use around livestock.

KEY FEATURES INCLUDE:

- S.W.A.T. will provide over a months' worth of 24-hour fly protection, at 10-minute spray intervals
- Feature filled dispenser, with the option to dispense at 5, 10, 15 and 30 minute intervals, 24 hours per day or only during daylight or darkness
- The dispenser features a flashing LED to show when the unit has power, this means when the unit stops flashing the batteries need to be replaced
- A super compressed, smaller and more environmentally friendly can

PART NO	DESCRIPTION
549030	SWAT Auto Dispenser c/w SWAT Refill 300ml

S.W.A.T REFILL CAN FOR DISPENSER



This is a replacement aerosol for S.W.A.T, the automatic fly control system. Each refill will work for over a month, with a powerful, yet natural Pyrethrin based broad spectrum insecticide which is, non-toxic and safe to use around livestock.

PART NO	DESCRIPTION
549025	SWAT Refill Can 300ml For Dispenser

PIG SORTING PADDLE



The Pig Sorting Paddle is a humane alternative for livestock, pig sorting and handling. The brightly coloured paddle head is 31 cm long, 16 cm wide and acts as a visual guide.

KEY FEATURES INCLUDE:

- Small ball bearing pieces inside the head that make a sound to which animals respond.
- Brightly coloured, so acts a visual guide

PART NO	DESCRIPTION
77001	Pig Sorting Paddle

RITCHEY ANTICHEW LIQUID 500ML

The Ritchey Anti Chew is a quick drying astringent bitter spray, which dries on the area and reduces bleeding. Contains a purple colourant, to identify the pigs in which the Anti-chew has been applied.



KEY FEATURES INCLUDE:

- Curbs tail and flank biting in growing and finishing pigs
- Purple colourant makes the pig easily identifiable

PART NO	DESCRIPTION
M41401	Ritchey Antichew Liquid 500ml

TAR PLASTER



This quick drying spray forms a transparent waterproof film and is ideal for minor wounds, navels, cuts and scratches.

PART NO	DESCRIPTION
495098	Tar Plaster Aerosol 200ml

- * Sold as a single item

BRITISH HAND CRAFTED POCKET KNIVES



A range of quality handmade knives, crafted in England by apprentice trained Sheffield Cutlers, using methods passed down from generation to generation. Every part of the process, from blade blanking, hardening, tempering and grinding through to the final polishing is completed in the English factory, ensuring a consistent and high quality product. Natural hues in the wooden handles can make every knife seem unique.

LAMBSFOOT KNIFE SCOTTISH PATTERN



Handmade in England, a very handy small bladed knife, with brass linings on a large, easy to grip rosewood handle with nickel bolsters. The blade is made of C70 Carbon Steel and is 5cm (2") long. The closed length of the knife is 10cm. Carbon Steel blades keep their edge and sharpen well.

PART NO	DESCRIPTION
430067	Lambfoot Knife Scottish Pattern 5cm (2")

LAMBSFOOT KNIFE LARGE



Handmade in England, in a traditional style, with brass linings and nickel bolsters. The blade is 8cm (3") long with a beautifully finished wooden handle.

PART NO	DESCRIPTION
430080	Lambfoot Knife Large 10cm (4")

LAMBSFOOT KNIFE SMALL



Handmade in England, in a traditional style, with brass linings and nickel bolsters. The blade is 7cm (2.75") long with a beautifully finished wooden handle.

PART NO	DESCRIPTION
430079	Lambfoot Knife Small 9cm (3.5")

EVERYDAY POCKET KNIVES

Economical pocket knives, great for everyday use with good quality blades and fixings.

LAMBSFOOT KNIFE NARROW BLADE



Wooden handle with a straight 8cm (3") blade.

KEY FEATURES INCLUDE:

- Economical
- Good quality blade and fixings

PART NO	DESCRIPTION
430201	Lambfoot Knife Narrow Blade 8cm (3")

COX STRAIGHT BLADE PENKNIFE



Plastic Handle, 6cm (2.5") blade.

PART NO	DESCRIPTION
430205	Cox Straight Blade Penknife 6cm (2.5")

COX CURVED BLADE PENKNIFE



Plastic Handle, 6cm (2.5") blade.

PART NO	DESCRIPTION
430207	Cox Curved Blade Penknife 6cm (2.5")

ANIMAL HANDLING & WEIGHING EQUIPMENT

We distribute both Prattley and Alligator mobile sheep and cattle handling equipment as well as Tru-Test Weighing equipment. In addition to mobile sheep and cattle yards our range also includes; three point linkage sheep yards, weigh crates, electric and auto drafters, handling race, weigh cells, load bars and weighing platforms. We also offer handling equipment for small holders.

CONTACT
THE ALLFLEX HUB
FOR MORE DETAILS
ON OUR PRATTLLEY
AND ALLIGATOR
WEIGHING AND
HANDLING SYSTEMS
HUB@ALLFLEX.CO.UK



Prattley

Whether it is a relocatable stock yard, a cattle crush or an automated drafting system Prattley has the answer.

All Prattley equipment is fully portable and can be taken to any part of the farm even in the roughest terrain allowing for on the spot stock management solutions.

Prattley Systems are also extremely tough and durable, for example the alloy used in production is heat treated which provides high tensile strength whilst still being light. This means that each system can be easily set up and operated by just one user.

FIND OUT MORE

For the full range of Prattley Animal Handling Equipment contact your local account Manager, visit www.ritchey.co.uk or email the Hub hub@allflex.co.uk



TAKE THE PENS TO THE SHEEP NOT THE SHEEP TO THE PENS

Moving sheep around the farm takes up valuable time, causes stress and poses an increased risk to biosecurity. Therefore keeping your livestock movements to an absolute minimum helps to improve productivity whilst reducing the risk of biosecurity.

Alligator has been helping British farmers for over four decades and has changed very little since it was originally launched. It's a system which is lightweight, modular, engineered for strength and above all else is designed to be fully mobile.

Because every farm is different farming large areas and split sites can be challenging. Alligator is easy to use and straight forward to set up what is more it can be erected by just one person. The system can handle up to 400 ewes and with the option to add up to three shearing pens it can be set up and operational in less than 30 minutes.

Alligator allows you to **Treat / Draft / Weigh / Footbath / Shear / Crutch and Handle**

The system is:

- ✓ **Mobile** Alligator is lightweight, robust and can be erected in minutes
- ✓ **Durable** Alligator boasts high grade raw materials, design and build quality
- ✓ **Versatile** Totally modular and reconfigured in minutes, Alligator is an ideal solution to remote sheep handling
- ✓ **Easy to use** Super strong and lightweight Alligator is exceptionally easy to use, even single-handedly
- ✓ **Cost effective** Built to last saving time and money

The basic Alligator system is the most popular allowing up to 200 sheep to be handled at any one time. Many choose to add a diamond forcing pen to this, running smaller numbers from the large gathering pen into this and then from there into the race for treating or drafting as required. Additional components can be added as and when required such as dosing race, weighing systems, shearing pen and footbath etc.

Find out more:

Our technical support team can help advise you on the right Alligator package to suit you. What's more we can offer a free farm demonstration to help you in deciding which system is best suited to your needs. For example, a basic system may be your starting point with component additions - all held in stock and ready for immediate delivery. Contact your local account manager or email us at hub@allflex.co.uk



ELECTRONIC WEIGHING

TRU-TEST®

WHY WEIGH?

Weighing replaces guess-work with facts, making it a powerful performance measurement tool for farmers.

The information provided by Tru-Test weigh scale indicators assists in decision making about feed rations, health treatments and condition scoring. On its own, weighing is a valuable addition to your operation. However, with the emergence of electronic identification (EiD), an opportunity exists to maximise the potential of the data collected.



WHAT IS EiD?

Electronic Identification (EiD) allows you to track individuals through their entire life cycle via an electronic ear tag (or other EiD device such as a bolus).

Many governments have implemented (or are planning to implement) mandatory EiD schemes to help protect the industry from the impact of biosecurity threats. From a farmer's point of view, tagging animals with an electronic (EiD) tag can add value by allowing better management of individual performance.

With EiD, an animal's history can be tracked allowing identification of good and bad performers. Streamlining your herd/flock allows inputs to be assessed and profits maximised.

For details of EiD tags see the Livestock Identification section.

WHAT ARE THE BENEFITS OF WEIGHING TO YOUR FARMING OPERATION?

Avoid Penalties! Ensure livestock are sold at the optimum weight grade

Closely managing the weight gain of individual animals or groups can help you tailor your feed and treatment inputs to ensure maximum profitability, reduce the risk of penalties and saves you money!

Minimise treatment waste

Drench/pour-on dosage can be calculated more accurately based on the requirements of the individual animal. Incorrect dosing is the number one cause of drench resistance amongst livestock, so while minimising treatment waste, you're also maximising the health of your herd.

Monitor reproductive performance

Closely monitoring the weight gain of offspring can help make breeding decisions to maximise reproductive performance.

In addition, ensuring optimum mating weights can improve breeding performance. Making sure target calving weights are reached will improve calving outcome.

Identify sick or lame animals earlier

Regular weighing of individual animals means weight loss and other issues that manifest over time are identified sooner, allowing prevention rather than cure. Early detection helps when managing input costs such as feed, health treatment and culling.

Manage your time spent in the yard

Your time is a valuable input in your farming operation. With the use of Tru-Test weigh systems your weigh sessions are completed faster and with less labour.

WEIGH INDICATOR SELECTION CHART



EziWeigh 5i



EziWeigh 7i



ID5000



XR5000

The EziWeigh5i Weigh Indicator is simple to use, easy to read, accurate and robust and is just what you need when collecting and viewing weights is your priority

Efficiently record and track individual animal weight gain. The EziWeigh7i provides the facts to help you manage your operation for increased profitability

Quickly capture accurate animal weight data to maximise profitability. Use it to identify poor performers and catch health issues early

Capture a wide range of animal information to make the right decisions. Quickly identify poor performing animals, optimise breeding programmes and catch animal health issues early

Part code	61513	61515	61310	61309
Accessories included	Mains Power adaptor, USB cable and mounting bracket	Mains power adaptor, USB cable, USB flash drive, mounting bracket, Data Link computer software (supplied on USB drive)	Hard plastic carry case, mains power adaptor, USB cable, USB flash drive, adjustable mounting bracket, EziLink computer software (supplied on USB drive)	Hard plastic carry case, mains power adaptor, USB cable, USB flash drive, adjustable mounting bracket, EziLink computer software (supplied on USB flash)
Size	205 mm (w) x 220 mm (h) x 74 mm (d)	205 mm (w) x 220 mm (h) x 74 mm (d)	260 mm (w) x 225 mm (h) x 65mm (d)	260 mm (w) x 225 mm (h) x 65 mm (d)
Weight (includes packaging)	2.2kg	2.2kg	4.6kg	4.6kg
Rechargeable Battery	✓	✓	✓	✓
Display	Large 35 mm display	Large 35 mm display	7" Sunlight viewable transfective LCD	7" Sunlight viewable transfective LCD
Backlit Screen	✓	✓	✓	✓
Battery Life and charge time	20 hours / 3 hours' charge time	15 hours / 3 hours' charge time	13 hours with backlight off / 5 hours' charge time	13 hours with backlight off / 5 hours' charge time
Computer Link - Desktop - Windows only - Mobile device (Andriod or Apple IOS)	✗	✓	✓	✓
Plug and Play	✓	✓	✓ predefined setups	✓ predefined setups
Recording & Memory	None, Live weight display only	Stores up to 20,000 records in 100 weighing sessions	Stores up to 250k records in 250 weighing sessions	Stores up to 1 million records in 1000 weighing sessions
EID Compatibility	✗	✓	✓	✓
Bluetooth Connectivity	✗	✓	✓	✓
Comment Field	✗	Yes - can record notes for each individual animal	Session Statistics / Animal Weight Summary / history table	Session Statistics / Animal Weight Summary / history table
Super Damp III Technology	✓	✓	✓	✓
Daily Live Weight Gain Display	✗	✓	✓	✓
Group Statistics	✗	✓	✓	✓
Individual Statistics	✗	✓	✓	✓
Auto Draft Compatible	✗	✓	✓	✓
App – Apple and Android	✗	✓	✓	✓
Withdrawal Period Alert	✗	✗	✗	✓

All indicators are IP67 rated and come with a 24-month warranty against manufactural defect

ELECTRONIC WEIGHING



ABOUT IP67

IP Ratings: (Ingress Protection or International Protection) - The Ingress Protection rating system is a classification system showing the degrees of protection from solid objects and liquids. The first number refers to the protection against solid objects, normally dust in relation to testing Pocket PCs. If the first number is 0, there is no protection provided. A number 5 refers to limited protection against dust. The number 6 is for total protection against dust.

The second number of the IP rating system refers to protection against liquids. A '0' indicates no protection, while a 7 refers to protection against immersion between 15 centimetres and one meter.

MP SERIES LOAD BARS



Tru-Test load bars are recognised for their rugged durability and performance in harsh environments. The internal cells are made from aircraft grade aluminium and are fully sealed to keep out moisture. Heavy duty TPR/PVS cables are a standard feature, which minimises the likelihood of cable damage.

Heavy Duty load bars or the Extra Heavy Duty load bars should be used for cattle crushes / crates.

The MP load bars are versatile and can be used for all animal types they can be used under cattle platforms, sheep crates and wool presses.

In-House Servicing and Repair available
Contact us to find out more

LOAD BAR SELECTION CHART	MP600	MP800	MP1010	HD800	HD1010
Part code	60001	60101	60201	60102	60202
Total capacity	2,000 kg	3,000 kg	3,000 kg	3,000 kg	2 x HD load bars = 6,600 lbs (3,000 kg)
Total Nominal length	600 mm	800 mm	1010 mm	800 mm	1,010 mm
Weight (pair)	13.5 kg	16.5 kg	22 kg	22 kg	39 kg
Cable length	5m	5m	5m	5m	5m
Resolution: Fine Mode	0.1 kg (0.2 lb) when used with Tru-Test weigh scale indicators			Base resolution with Tru-Test indicators: 2 x HD load bars = 0.5 kg	Base resolution with Tru-Test indicators: 2 x HD load bars = 0.5 kg
Resolution: Normal	0.5 kg (1lb)				
Accuracy	+1% or two resolutions, whichever is greater			+1% or 2 resolutions, whichever is greater	+1% or 2 resolutions, whichever is greater

MP MULTI-PURPOSE LOAD BARS

KEY FEATURES (MP600 MP800 MP1010)

- **Versatile:** Can be temporarily installed and when used with a lightweight Tru-Test aluminium platform they can be easily transported between different weighing sites.
- **Compatible:** They can be used with any weigh scale indicator in the Tru-Test range without any need for calibration.
- **Corrosion resistant:** Zinc plated gold passivated load bar housing for additional corrosion resistance.
- **Sealed load cell:** Protects the load cell against moisture damage allowing the load bars to remain outside.
- **Easy set up:** Simply requires a firm level surface and a rigid crate or platform.



HD HEAVY-DUTY LOAD BARS

With a 6,600-lb capacity, the HD load bars are especially suited for use under manual squeeze crushes. Their industrial strength construction makes them the load bar of choice for cattle operations where there are severe loadings and high usage rates.

KEY FEATURES (HD800 AND HD1010)

- **Permanent/static installations:** Designed for installation under manual squeeze crushes
- **Compatible:** Can be used with any weigh scale indicator in the Tru-Test range without calibration
- **Corrosion resistant:** Galvanized steel housing protects the bar in the tough agricultural environment
- **Sealed load cell:** Protects the load cell against moisture damage allowing the load bars to remain outside under the crush
- **Heavy duty:** Designed for high volume, large weight use within a crush environment
- **Strong and durable construction:** Ensuring long term reliability



SUSPENSION CELL

Versatile, rugged, ideal for mechanical scales conversions. The Tru-Test Suspension Cell allows you to weigh stock accurately and quickly due to the simple, robust yet highly sophisticated EziWeigh indicator allied to the suspension cell.

KEY FEATURES

- 250 kg capacity
- 2 mm threaded mounting holes



PART NO	DESCRIPTION
60612	Suspension Cell

All indicators are IP67 rated and come with a 24-month warranty against manufacturing defect

EZIWEIGH5i

The EziWeigh5i Weigh Indicator is simple to use, easy to read, accurate and robust and is just what you need when collecting and viewing weights is your priority.

KEY FEATURES

- **Highly visible display:** Large 35 mm display with a super bright backlight making it viewable in all conditions, from brightest daylight to darkness.
- **Plug and play load cell connection:** No calibration is required as the weigh scale automatically recognises the Tru-Test load sensing device connected.
- **Superdamp™ III technology:** Capture accurate live weights more quickly than any other weigh scale indicator on the market. The EziWeigh5i has improved weighing electronics to lock onto lively animals even faster.
- **Auto Zero:** Option to ensure any accumulation of animal dirt on the platform or crate is zeroed out.
- **Long life and fast charging battery:** Up to 20 hours of run time and only 3 hours to recharge. Comes with mains power adaptor and USB.
- **Durable case:** IP67 waterproof case with double walled design for extra protection.



PART NO	DESCRIPTION
61513	EziWeigh5i

EZIWEIGH7i



WEIGH, TRACK, ACT.

The EziWeigh7i allows you to efficiently record and track individual weight gain. Providing hard facts to help you manage your operation with greater certainty and profitability. With the capability to record information against Electronic ID (EiD) or Visual tags, you can easily capture and analyse individual animal performance and history.

Replacing guesswork with facts, the EziWeigh7i displays live weight gain since the previous session, showing animal performance right there in the yard. Using its built-in Bluetooth® technology, the EziWeigh7i wirelessly connects to Tru-Test's EiD readers, automatically transferring scanned ID's straight into the weigh scale indicator, saving time and improving accuracy. Alternatively, the animal's VID can be entered using the large and durable keypad.

KEY FEATURES

- **Wireless or tethered EiD reader connectivity:** Choose between wireless Class 1 Bluetooth® communication with EiD readers up to 100 m away, or connect via cable.
- **ID capable:** Records information against Electronic (EiD) and Visual tags. VID numbers can be automatically cross-referenced with their EiD match, saving you time and reducing data entry.
- **Live weight gain:** Calculates and displays an animal's daily weight gain since the previous session.
- **Capture traits or treatments:** Record notes for each individual animal.
- **Memory and storage capacity:** Stores up to 100 weighing sessions and 20,000 records.
- **Easily transfer recorded information:** No matter where you are you can quickly transfer information using your mobile device (Android or Apple iOS) or computer (Microsoft Windows).
- **Superdamp™ III technology:** Capture accurate live weights more quickly than any other weigh scale indicator on the market. The EziWeigh7i has improved weighing electronics to lock onto lively animals even faster.
- **Statistics:** Summary statistics from the current and previous weigh sessions can be recalled right there in the yard.
- **Long Battery Life:** Over 15 hours of run time.



- Auto incrementing and auto prefixing.
- Bright red stable light.
- IP67 rated, 100% water and dust proof.
- Android and Apple Data Link Apps compatibility.

**In-House
Servicing and
Repair available**
Contact us to
find out more

PART NO	DESCRIPTION
61515	EziWeigh7i

THE 5000 SERIES

The top of the range 5000 series featuring the ID5000 and the XR5000 weigh scale indicators are the ultimate in powerful business management data tools allowing you to capture a wide range of key animal information. Both allow you to quickly identify poor performing animals, optimise breeding programmes and catch animal health issues before they become problematic.

KEY FEATURES

- Large 7" display with unique transfective screen technology
- Fast, intuitive navigation
- Huge amount of data storage
- Simple step-by-step setup
- Ability to save favourite setups
- Water, shock and dust proof, IP67 rated
- Android and Apple data Link apps

PART NO	DESCRIPTION
61310	The ID5000
61309	The XR5000



AT A GLANCE THE ID5000 AND XR5000 WEIGH SCALE INDICATORS

	ID5000	XR5000
How many records can I store?	250,000	1 million
How many weighing sessions can I record?	250	1,000
How many pieces of data can I keep on each animal?	3 (Fixed)	100 (Custom)
How many tag identifications can I assign to an animal?	3	4
Can I customise the information that I record?	✗	✓
Can I view animal history statistics in the yard?	✓	✓
Can I capture an animals weight gain?	✓	✓
Can I draft by weight gain or single criteria?	✓	✓
Can the indicator control an AutoDrafter?	✓	✓
How many ways can I draft?	5	10

THE ID5000



THE XR5000



EID XRS2 STICK READER



DATA COLLECTION MADE EASY

The Tru-Test XRS2 Stick Reader is the next generation of readers and offers truly outstanding data recording capacity allowing you to scan your animals rapidly. The XRS2 provides you with key information to be able to make informed and accurate management decisions on individual animal performance and back up critical health, breeding and culling decision in real time.

The XRS2 Stick Reader is suitable for cattle, sheep and deer, and features a large highly readable screen and alphanumeric keypad for capturing important information from the pasture or pens.

Coupled with Bluetooth® wireless technology it allows you to connect your reader to your Weigh Scale Indicator and means you can easily share information using your smartphone.

PART NO	DESCRIPTION
829955	EID XRS2 Stick Reader

KEY FEATURES

- Record up to 1 MILLION records
- Alerts for pre-selected animals
- Large sunlight viewable 2.7" colour LCD screen making it easy to read
- Alphanumeric keypad for flexible data entry
- Fast tag reads with feedback through vibrating handle
- Customisable alerts such as cull
- Add and edit session names to keep track of each session
- Capture animal information such as condition scores in up to 15 custom data fields
- Easily download recorded information using your Windows computer or Apple iOS/Android smartphone
- Up to 19 hours read operation and 9.5 hours of continuous read
- Long reach

THE ERS HANDHELD READER



DATA COLLECTION MADE EASY

The Tru-Test ERS Handheld Reader has been designed to eliminate the hassle of reading and recording livestock tags manually, Electronic Identification (EiD) ear tags can be easily logged by simply waving the reader by the animal's ear.

It is compact and lightweight with an ergonomic rugged waterproof and dustproof casing and what's more its IP67 rated which makes it perfect for tough agricultural environments. More than capable as a standalone EiD tag reader and with enough memory to store up to 4000 records with a battery life of 8 hours.

The ERS Handheld Reader can be easily linked up to a PC for downloading tag information, alternatively when paired with Tru-Test® weigh scale indicators via Bluetooth it becomes an integral part of your Complete Livestock Management Solution.

KEY FEATURES

- Compact, simple and easy to use
- 4,000 tag memory
- Fully rechargeable
- 8 hours' battery life
- Duplicate discard
- Bluetooth® enabled
- Connects with Tru-Test® weigh scale indicator and a PC



PART NO	DESCRIPTION
09790	The ERS Handheld Reader

MOBILE PRINTER WITH BLUETOOTH



KEY FEATURES

- Fully portable and compatible with the XRS2 Stick Reader
- Bluetooth for wireless communication
- Magnetic base
- Rechargeable
- Easy to load paper roll



PART NO	DESCRIPTION
09819	Mobile printer with Bluetooth

PANEL READERS

XRP2 PANEL READERS



The XRP2 is built around Tru-Test's state of the art digital signal processing that gives optimum read range while automatically reducing unwanted interference.

Its revolutionary fast auto tuning provides maximum performance under all conditions. The fixed panel reader environment isn't static. Small changes in the metal work or the simple action of opening a metal gate on a crush or crate can affect the reader's range and read speed performance. Tru-Test's EID Panel Readers maintain optimal performance through these changing conditions by retuning up to 18 times per second, ensuring consistent, accurate tag reads throughout your weigh session.

The XRP2 gives you the flexibility to record EID tags into a session file or send EIDs directly to your weigh scale indicator to record individual animal weights and other information. It also comes with built-in Bluetooth® wireless technology to easily connect to your weigh scale indicator or smartphone.

KEY FEATURES

- Reads all ISO HDX and FDX-B tags
- Option to record tags or send to your weigh scale indicator
- Wirelessly connect with your weigh scale indicator or smartphone using Bluetooth® technology
- Easily download recorded information using your PC or Android smartphone
- Plug and play setup, with unique auto tuning function
- Large and small antenna options to suit your setup
- Rugged design for the tough agricultural environment



PART NO	DESCRIPTION
09794	XRP2 Panel Reader - Sheep
09793	XRP2 Panel Reader - Cattle



LARGE ANTENNA

- Rugged design and easy to install
- Large antenna
- Read range single antenna up to 39" (1 m)
- Antenna size: 24x34 inches / 58.42x86.36cm (Width x Height)
- Weight: 18 lb / 8kg

SMALL ANTENNA

- Rugged design and easy to install
- Small antenna
- Read range single antenna up to 30" (0.75m)
- Antenna size: 16x17 inches / 40.64x43.18cm (Width x Height)
- Weight: 5 lb / 2kg

WORLD CLASS PRECISION MADE CLIPPING AND SHEARING EQUIPMENT

SWISS ENGINEERING AT ITS BEST

The Heiniger range of clippers, hand pieces and blades combine swiss engineering and outstanding blade technology for a perfect finish every time. Providing exceptional performance including reliability, ease-of-use, clip finish and value for money.

All Heiniger clippers, hand pieces and blades come with a two-year warranty as standard and what is more with our Customer After Sales Service department we can service and repair all Heiniger clippers as well as provide a blade sharpening service all under one roof.

Key features:

- Powerful and fast
- Lightweight ergonomic design
- Comfortable
- Reliable and quiet
- 2 Year warranty

**IN-HOUSE
SERVICING
AND REPAIR
AVAILABLE**

Contact us to
find out more

CLIPPING AND GROOMING AT A GLANCE: EQUINE / CATTLE / SHEEP

MAINS POWERED UNITS

HEINIGER PROGRESS CLIPPER



The Heiniger Progress Clipper has been designed with a slim grip is lightweight and quiet in operation.

The Clipper is the best solution to clip nervous horses and uses an innovative permanent magnet motor which provides a powerful but quieter clip.

PART NO	DESCRIPTION
538038	Heiniger Progress Clipper - Clipper fitted with 31F - 15 Fine Blade Set. Includes a spare blade set.

HEINIGER PROGRESS STYLE



The Progress Style Clipper is a trendy design clipper. Equipped with the latest electronic technology, the Progress Style horse clipper is now even more powerful and quieter than ever.

Weighing in at only 930 grams, it has an impressive slim ergonomic handle making it a popular choice for people whom wish to clip their own few horses at a time that suits them.

PART NO	DESCRIPTION
538040	Heiniger Progress Style - Clipper fitted with 31 - 15 Standard Blade Set.

HEINIGER XPERIANCE CLIPPER



The Heiniger Xperience Clipper combines power, performance and an ergonomic lightweight design offering a unique aspect in the innovative, stable connection between the clipping head and the motor housing.

PART NO	DESCRIPTION
538044	Heiniger Xperience Clipper - Clipper supplied with a 21 - 23 Coarse Blade Set. Includes a spare blade set.


HEINIGER XPERT CLIPPER



The Xpert sheep shearing machine is a real power package. Power and technology are combined at the highest level in a handy, lightweight machine. The unique and innovative, strong connection between the shearing head and the fibreglass-reinforced housing ensures low vibration shearing, with an ergonomic handle design.

PART NO	DESCRIPTION
538050	Sheep Clipper - Clipper supplied with Ovina Premium Comb and Jet Cutter. Includes a spare Comb and Cutter.


HEINIGER CATTLE CLIPPER 12 VOLT



A versatile clipper for use out in the field and in remote locations, combining strength with comfortable handling and ease of use. The clipper can operate for 4-5 hours from a car battery and it can also be converted to a sheep clipper with a sheep head.

PART NO	DESCRIPTION
538177	Cattle 12 Volt Clipper - Clipper supplied with a 21 - 23 Coarse Blade Set. Includes a spare blade set.


HEINIGER SHEEP CLIPPER 12 VOLT



A versatile clipper for use out in the field and in remote locations. Glass fibre reinforced polyamide casing combines strength with comfortable handling and ease of use. The clipper can operate for 4-5 hours from a car battery and is ideal for dagging and crutching.

PART NO	DESCRIPTION
538054	Sheep S12 Volt Clipper - Clipper supplied with a Ovina Premium Comb and Jet Cutter. Includes a spare comb and cutter.

HEINIGER SAPHIR CORD



The Saphir Cord is based on the established Saphir Clipper. The Saphir Cord is three metre mains corded unit, which gives the clipper 3200 double blade strokes per minute, allowing it to work easily through all types of coats. The unit is supplied with a 1.8 mm No 10 blade.

PART NO	DESCRIPTION
539410	Saphir Cord

HEINIGER HANDY CLIPPER *



The Handy is a quiet lightweight clipper weighing only 1.3 Kg and is therefore ideal for clipping nervous animals. The slim body is perfect for small hands and is easy to use in tight areas. The powerful 120 watt cooled motor delivers a high cutting speed of 3200 rpm for speedy clipping.

PART NO	DESCRIPTION
538018	Heiniger Handy Clipper - Clipper fitted with a 31F - 15 Fine Blade Set. A spare 31F- 15 Fine Blade Set is included for continuous clipping.

CORDLESS UNITS

HEINIGER XPLORER CORDLESS CLIPPER



BEST SELLER

A cordless clipper which uses Lithium-Ion battery technology providing up to two-hour cordless clipping from a single charge. No disturbing cable increases the flexibility and portability of the clipping experience. This X-Series clipper is very powerful, ergonomic and quiet. Through the unique double tooth gear, the power is transmitted optimally from the motor to the clipping head. The practical LED display on the charger unit indicates the current charge level of the battery.

Incl. carry case, cutter, screwdriver, oil, brush.

PART NO	DESCRIPTION
538046	Heiniger Xplorer Cordless Clipper - Clipper supplied with 31F - 15 Fine Blade Set. Includes a spare blade set.

ALSO AVAILABLE
Xplorer Cordless Clipper c/w Extra Battery (part no XB538046)

HEINIGER SAPHIR CORDLESS CLIPPER




Up to 90 minutes clipping time per battery with high performance lithium-ionic technology. It is supplied with an additional battery pack and less than 60 minutes charging time for continuous clipping.

The perfect clipping machine for trimming cattle and horses or cow-styling. The detachable Saphir blade system means it is compatible with other clippers such as Andis, Moser, Oster and Wahl.

PART NO	DESCRIPTION
539210	Heiniger Saphir Cordless Clipper - c/w No 10 Blade.

HEINIGER SAPHIR STYLE CORDLESS CLIPPER



The perfect clipping machine for trimming horses or for horse styling. The Saphir Style has low noise but delivers high performance from its 35 W permanent magnet motor generated from its Lithium-Ion battery for continuous clipping. The Style is supplied with a high-quality charger and an additional charging dock for a second battery and charging indicator. The style also comes with a large selection of interchangeable clipping blades which are compatible with Andis, Moser, Oster and Wahl Clippers.

PART NO	DESCRIPTION
539310	Heiniger Saphir Style with no 10 blade

HEINIGER MIDI TRIMMER




IDEAL FOR DOG GROOMERS

This trimmer is used for the part trimming of a small animal and the contours of the head, legs, belly and tail. The adjustable cutting head allows variable cutting lengths between 0.5 - 2.5 mm. The package also includes five Snap-on combs in cutting lengths of 4, 8, 12, 16 and 20 mm. Furthermore, the clipper comes with an adapter allowing anytime trimming from 230 V mains. Light with low-noise and a precise cut, the Heiniger Style Midi Trimmer glides through the coats with ease.

PART NO	DESCRIPTION
539500	Heiniger Midi Trimmer

HEINIGER MINI TRIMMER



**IDEAL
FOR DOG
GROOMERS**

The Mini Trimmer trims quietly and with precision around the paws, eyes and tail of a small animal such as a dog. The adjustable Snap-On comb allows a cut length of 3 to 6 mm. This trimmer is ideal for the final styling of a dog. After trimming the trimmer head can easily be removed and cleaned.

PART NO	DESCRIPTION
539600	Heiniger Mini Trimmer

EQUINE / CATTLE AND SMALL ANIMAL BLADES

Heiniger blades are manufactured to the highest quality and specification to give optimum performance and maximum usage Swiss Engineering blade technology at its best. Blades sold at set.

STANDARD BLADE SET



A General purpose blade set for horses and cattle.

- Coat length after clipping 2-4mm
- Bottom blade 31 tooth - Top blade 15 tooth

PART NO	DESCRIPTION
538071	Standard Blade Set 31-15

COARSE BLADE SET



This blade is used particularly for beef and dairy cattle.

- Coat length after clipping 2-4mm
- Bottom blade 21 tooth - Top blade 23 tooth

PART NO	DESCRIPTION
538097	Coarse Blade Set

FINE BLADE SET



A general purpose blade set for horses and cattle.

- Coat length after clipping 1-2 mm
- Bottom blade 31F – Top blade 15 tooth

PART NO	DESCRIPTION
538102	Fine Blade Set

THOROUGHbred BLADE SET



An alternative to the standard with extra top blade teeth.

- Bottom blade 31 tooth – Top blade 23 tooth

PART NO	DESCRIPTION
538188	Thoroughbred Blade Set

SURGICAL BLADE SET



A surgical blade used for very close clipping such a veterinary work.

- Coat length after clipping 0-1 mm
- Bottom Blade 53 - Top blade 23 tooth

PART NO	DESCRIPTION
538117	Surgical Blade Set

SAPHIR BLADES

(Not Illustrated)

PART NO	DESCRIPTION
539107	Saphir Blade No.7F 3.2mm for Dogs / Horses
539110	Saphir Blade No.10 1.8mm for Dogs / Horses
539115	Saphir Blade No 15 1.5mm for Dogs / Cats
539140	Saphir Blade No 40 0.6mm surgical blade

HEINIGER SHEARING EQUIPMENT



HEINIGER ICON CYCLONE HAND PIECE



The Heiniger Icon Cyclone hand piece is designed using the latest state of the art technology resulting in a smoother cooler running hand piece, providing a crisper cutting action whilst supporting lower maintenance and an extended working life.

PART NO	DESCRIPTION
541202	Heiniger Icon Hand Piece (Worm Drive)
541200	Heiniger Icon Hand Piece (Pin Drive)

PLEASE NOTE: Ferrule sold separately.

EVO 3 SHEARING MACHINE

The first machine to be manufactured with an electronic safety switch which will stop the machine immediately if a lock up occurs. The EVO has three speed settings from 2700, 3200 and 3500 rpm allowing it to cater for all shearing conditions.

The self-tensioning belt means less maintenance and less wear and tear. Machines are available with solid or flexi drives with either worm or pin.

PART NO	DESCRIPTION
541614	EVO Machine complete with Flexi Worm Drive
541618	EVO Machine complete with Solid Worm Drive
541620	EVO Machine complete with Solid Pin Drive

EVO 3 MACHINE ONLY

KEY FEATURES

- Dimensions (WxHxL): 225 x 300 x 290mm
- Weight 7.8kg
- Noise emission 70dB

PART NO	DESCRIPTION
541630	EVO 3 Machine only

EVO ONE SHEARING MACHINE

Weighing only 6.5kg, the ONE is a single speed machine which incorporates many of the safety features found in the EVO. With a running speed of 2950rpm, the ONE can only be operated with flexi drive and is an ideal machine for farmers.

ONE MACHINE ONLY

KEY FEATURES

- Dimensions (WxHxL): 243 x 317 x 209mm
- Weight 6.5kg
- Noise emission 70dB

PART NO	DESCRIPTION
541650	ONE machine only

EVO/ONE FLEXI PIN DRIVE ONLY

PART NO	DESCRIPTION
541634	EVO/ONE Flexi Pin Drive only 200cm

EVO/ONE FLEXI WORM DRIVE ONLY

PART NO	DESCRIPTION
541632	EVO/ONE Flexi Worm Drive only 200cm

SOLID WORM DRIVE

PART NO	DESCRIPTION
541512	Solid Worm Drive

SOLID PIN DRIVE

PART NO	DESCRIPTION
541514	Solid Pin Drive

541650

FERRULE RIGID DOWN TUBE

PART NO	DESCRIPTION
542078	Ferrule Rigid Down Tube

FERRULE FLEXI DOWN TUBE

PART NO	DESCRIPTION
542080	Ferrule Flexi Down Tube

TALLY COUNTER

0 - 9999 push button operation, chrome plated with instant reset.



PART NO	DESCRIPTION
422010	Pocket Tally Counter



541630



Hand pieces sold separately

HEINIGER SHEARING EQUIPMENT - COMBS AND CUTTERS

A large selection of Combs and Cutters are available for shearing for open woolled sheep through to tough shearing crossbreds. Contact your local Cox & Ritchey Account Manager or our customer services department to find out more.

FARMER'S PACKS

For many, the ideal way to purchase cutters and combs. Each pack contains 2 cutters and 1 comb in a display pack.

STANDARD FARMER'S PK.



2 Jet cutters and 1 Ovina comb.

PART NO	DESCRIPTION
541080	Standard Farmer's Pack

WIDE FARMER'S PK.



2 Jet cutters & 1 Pro Legend comb.

PART NO	DESCRIPTION
541082	Wide Farmer's Pack



Trafalgar Hand Shears are a traditional Double Bow pattern made from high carbon steel.

KEY FEATURES

- Hand forged, high quality blades, hardened and tempered
- Hollow ground and hand set by blade specialists, to give optimum performance
- Attractively packaged for retail display

SMALL DAGGING SHEARS



Strong cut, handy pocket size, ideal for smaller hands.

PART NO	DESCRIPTION
418140	Small Dagging Shears 90mm (3.5")

STANDARD DAGGING SHEARS



Strong cut, traditional size, straight handles.

PART NO	DESCRIPTION
418142	Standard Dagging Shears 135mm (5")

BENT HANDLE DAGGING SHEARS



Strong cut, traditional size, bent handles.

PART NO	DESCRIPTION
418144	Bent Handle Dagging Shears 135mm (5")

STANDARD SHEARING SHEARS



Firm squeeze, long blades, good wool penetration, straight handles.

PART NO	DESCRIPTION
418146	Standard Shearing Shears 165mm (6.5")

BENT HANDLE SHEARING SHEARS



Firm squeeze, long blades, good wool penetration, bent handles.

PART NO	DESCRIPTION
418148	Bent Handle Shearing Shears 165mm (6.5")

CLASSIC SHEEP SHEAR SHARPENER

Handy pocket sized sharpener.

PART NO	DESCRIPTION
418150	Classic Sheep Shear Sharpener



SIMPLE TO USE AND EFFECTIVE

HEAVY DUTY COTTON ROPE CATTLE HALTER



A chunky, heavy weight cotton rope, making it comfortable to handle. Length 3.46m, diameter 2cm

KEY FEATURES INCLUDE:

- Chunky and heavy weight
- Heavy weight cotton for comfort

PART NO	DESCRIPTION
406223	Heavy Duty Cotton Rope Cattle Halter

CIRCUITEER II CATTLE BLOW DRYER

The Circuiteer Cattle Blow Dryer is portable and lightweight and can be suspended, carried, or set on the ground, making it the perfect accessory to have around the farm. The blow dryer has a twin high-speed 2.1hp motors, and can be used with only one motor for sensitive areas.

KEY FEATURES INCLUDE:

- Supplied with 20-ft electric cord and a 15-ft air hose with blower tip and filter
- Equipped with 240V utility outlet
- Circuiteer cattle blow drier can power clippers and lights
- The blower uses no heating element
- Portable and lightweight



PART NO	DESCRIPTION
12000	Circuiteer II Cattle Blow Dryer

GROOMING SCISSORS



Stainless Steel

PART NO	DESCRIPTION
511134	Dressing Scissors 150mm (6")
511146	Dressing Scissors 180mm (7")

CLIPPING SCISSORS



Curved blade and cranked handles
180mm (7")

PART NO	DESCRIPTION
511018	Clipping Scissors

MARKING SCISSORS



PART NO	DESCRIPTION
511031	Straight blades and cranked handles 190mm (7.5")
511043	Curved blades and cranked handles 190mm (7.5")

PERTH COMB



The Perth Comb is perfect for everyday use on cattle as well as during the show season. The comb can be used to take dirt out of an animal's coat, as well as dealing with stubborn areas in longer hair.

KEY FEATURES INCLUDE:

- Pin tooth 125mm (5")

PART NO	DESCRIPTION
517070	Perth Comb 230mm (9")
517069	Perth Comb 125mm (5")

CURRY COMBS



PART NO	DESCRIPTION
517033	Nickel plated teeth both sides
517021	Nickel plated teeth one side
518025	Norton Curry Comb

CUSTOMER AFTER SALES SERVICE

SERVICING TEAM

WE OFFER IN-HOUSE SERVICING AND REPAIR (INCLUDING AUTHORISED WARRANTY WORK) ON:

- ✓ Tru-Test Electronic Weighing Equipment and EiD Readers
- ✓ Heiniger, Oster and Andis Clippers and Shearing Equipment
- ✓ Gas and Electric Dehorners (Selected models only - contact us to find out more)
- ✓ Hotshot Pulsar Coaxers
- ✓ Phillips Drenchers and Vaccinators

We can also co-ordinate the repair of Prattley and Alligator Animal Handling Equipment

Turnaround timescales can vary; however, we will notify you during your initial enquiry.



SHARPENING SERVICE

We also provide a complete sharpening service for a wide range of clipper and shearing blades used by companion animal, agricultural, equine and veterinarian sectors and sharpen and realign scissors (Straight edge and Blenders) to the highest standards.

We aim to:

- ✓ Sharpen and dispatch blades and scissors within 2 to 3 working days upon receipt*
- ✓ Service/repair and dispatch all other items within 3 to 5 working days upon receipt*



CUSTOMER SERVICE

We strive to provide a first-class service and carry out our work as quickly and as cost effectively as possible for all our customers. Our friendly, experienced and knowledgeable team are always available to offer help and advice when it comes to keeping your equipment in the best possible working condition.

***Please note:** Subject to requested quote being accepted by the customer. The turnaround timescales can vary dependent on the volume of work at the time of receipt of enquiry.




Livestock Identification



ALLFLEX TAG RANGE

PRODUCTS BELOW SHOWN AS ACTUAL SIZE



COLOUR GUIDE



Colours may vary, please use as a guide only. Contact Customer Services to order free samples.



AN EXPANDING ROLE FOR THE HUMBLE EAR TAG

The collection, monitoring and archiving of animal data provides intelligence on key management parameters, such as performance, efficiency, welfare and genetic progress. Use of intelligence in this way is increasingly a key ingredient in efficient, profitable and sustainable production. At one time, just a simple method of animal identification, the humble ear tag has evolved to play an increasingly significant role in this new area of livestock farming intelligence.

In-ear tag technology can help to improve livestock performance and profitability, especially when allied with integrated handling and weighing systems. Dairy farmers routinely use information provided by milking parlour software, such as milk yield and somatic cell counts, to make informed management decisions for breeding or marketing. In the sheep and beef sectors the use of EiD (Electronic Identification) tags linked to farm management software can offer the same opportunity to focus on those factors most likely to influence profitability such as daily liveweight gain and finishing weights.

Of course, information or intelligence on individual animals can take many forms and cattle ear tags can also play their part in the management and eradication of disease - for example Bovine Viral Diarrhoea (BVD). Tissue Sampling Tags (TST's) take a small biopsy of the ear tissue at the point of application and when analysed can provide information about the disease status of the animal. This can facilitate informed culling and breeding decisions by allowing Persistently Infected (PI) animals - those that carry and spread the disease - to be identified.

All of the livestock identification tags within the Allflex UK Group portfolio are available for supply through our network of Cox and Ritchey Trade Partners at competitive rates - please contact your Local Account Manager for pricing information and technical specifications.

"The humble ear tag has evolved to play an increasingly significant role in this new area of livestock farming intelligence."



BLUE HELIX TAGS HAVE NOW BEEN REPLACED WITH ALLFLEX ULTRA TAGS

CATTLE TAGS

ULTRA CATTLE TAGS

Given the legal requirement to identify cattle, livestock producers can gain significant commercial benefit from utilising in-ear tag technology to make informed management decisions. Using Tissue Sampling Tags can help to manage disease whilst EiD tags can help to make breeding decisions or market stock at the right weight and specification. Used by livestock producers across the world, Allflex Ultra cattle tags are manufactured from high quality plastic and incorporate design features to ensure free rotation whilst in the ear; ensuring reliability and readability.

- Officially approved and available as senior (large) and junior (medium) flag and button
- Flag tags have reinforced shoulders providing robustness
- 360° tag pin rotation allows tags to freely rotate in the ear
- Female tag is larger than male tag to allow for easier identification once applied
- Clear black laser printing for distance readability and longevity of the mark
- Available with management print and barcode options
- Available as HDX or FDX EiD or Tissue Sample Tags
- Apply with the B06 or B08 Applicator
- **Available colours:** Yellow, Blue, Green, Orange, Pink, Red, White and Purple

OFFICIAL PRIMARY TAGS



REPLACEMENT CATTLE TAGS

**ORDER BEFORE 3PM
FOR SAME DAY DESPATCH***


Express delivery options are available
- please contact Customer Services.

* Subject to BCMS approval.



OFFICIAL SECONDARY TAGS

SEE
PAGE 63
FOR TISSUE
SAMPLING
TAGS



MINIMUM ORDER QUANTITY 75


Manufactured from quality stainless steel, QWIKSILVER tags won't rust or corrode; ensuring readability and improved comfort as a result of lower rates of infection. The high performance applicator has been specially engineered and manufactured with high quality bearings for smooth application. Find the QWIKSILVER applicator on page 69.

CATTLE TAGS

ULTRA EiD CATTLE TAGS

- HDX (Half Duplex Tag) parlour tags are designed to work with all HDX parlour systems
- The technology optimises signal transmission and provides the greatest possible read distance
- For total control/traceability in the milking, feeding and identification of individual animals
- Choose between discreet EiD button tags, or EiD flag tags for greater visibility
- Apply with either the B06 or B08 applicator; remember to remove the black insert prior to fitting female EiD tags
- **Available colours:** Yellow and White

THE TAG BUCKET FILE WILL BE SENT VIA E-MAIL
Please supply a valid e-mail address



Applying EiD tags to calves at birth and using these in conjunction with an EiD reader and weigh indicator can significantly aid the monitoring and management of youngstock, allowing for informed management decisions.

Find the full range of Tru-Test EiD data recording and weighing equipment in our stock management solutions brochure



DID YOU KNOW?

RE-USEABLE EiD TAGS CAN BE USED AS NON-OFFICIAL MANAGEMENT TAGS

See page 64 to find out more.



CATTLE TAGS

TAG AND TEST FOR BVD FROM BIRTH

Tissue Sampling Tags (TST's) can be used as part of a management strategy to help eradicate Bovine Viral Diarrhoea (BVD). Doing so will help to increase overall productivity and profitability and could give a potential price premium through selling BVD free stock

- BVD is a highly contagious viral disease of cattle, most UK herds have already been exposed
- There can often be no obvious clinical signs of infection
- The disease can cause reduced growth rates and milk yield, abortion and secondary disease
- BVD is estimated to cost £46.50 per cow per year* and is one of the biggest issues facing the UK cattle industry


*Source: Yarnall and Thrusfield (2017) Vet Record doi: 10.1136/vr.104370

SEE PAGE 69 FOR THE TST APPLICATOR



An important tool towards a BVD free future...

- Officially approved cattle tags; use management TST's to test older cattle
- Apply with the B07MICRO TST Universal Applicator; applicators purchased prior to September 2017 may require an orange micro tube adapter (part code: 67001768)
- Take a sample of tissue as part of the normal tagging process
- The sample is tested in the laboratory to determine the BVD status of the individual animal
- Infected animals can be identified for culling
- BVD free animals can be recorded on your national database
- **Available colours:** Yellow and White



The official Ministry number is intrinsically linked to the tag, sampling needle and sample tube.

CATTLE MANAGEMENT TAGS

AVAILABLE EITHER PRINTED OR BLANK

Management tags can be used for a number of reasons and can assist in the identification, monitoring and ongoing management of selected animals.

Ideal for:

- ✓ Breeding decisions
- ✓ Pedigree management
- ✓ Disease management
- ✓ Cull selection
- ✓ Health status
- ✓ Temperament/behaviour characteristics



SENIOR
LARGE



JUNIOR
MEDIUM



EID BUTTON



REUSABLE OPEN
HEAD EID BUTTON

WITH LONG STEM MALE PIN
TO ENABLE EASY REMOVAL

4MM
LONGER
THAN A
STANDARD
MALE PIN



TST MANAGEMENT BUTTON

ALSO AVAILABLE AS SENIOR AND JUNIOR
MANAGEMENT IN WHITE OR YELLOW

Use TST management
button tags *applied with
a B07MICRO applicator for
testing older, pre-tagged
cattle.

TAG MARKER

An all weather long lasting, fade resistant
tag marker pen, supplied with a fine and
narrow nib.



SHEEP TAGS

EiD SLAUGHTER/BATCH TAGS



QWIKEID

MINIMUM ORDER 10

SLAUGHTER/
BATCH TAG
ONLY

A robust loop tag which incorporates a high quality transponder. Superior read range, high visibility laser marking and a precision-engineered applicator ensure this tag is reliable and easy to apply. Upon application the tag pin blunts.



KLIK TAG

MINIMUM ORDER 10

Easy-tear strips and a sharp, narrow pin allow for fast loading and tagging. The curved pin securely locks the tag on application to ensure retention.



RapID

MINIMUM ORDER 10

Speed up identification of your market ready stock with the RapID Semi-Automatic Tagging System – tag up to 300 lambs per hour! Supplied in strips of 20, the RapID can also be applied with the RapID Single Shot Applicator.

Please note this is a special order item, please contact us for lead times.

DID YOU KNOW?

We can print management information on slaughter/batch tags

SHEEP TAGS

EiD SLAUGHTER/BATCH TAGS



BUBBLEGUM

MINIMUM ORDER 10

Concealing the transponder in the tag pin allows for excellent weight distribution. Manufactured from ultra-flexible plastic, the tip of the tag pin has a cross cut point which helps tag application by cutting through ear tissue as opposed to punching through.



SNAPP

MINIMUM ORDER 10

Based on the original sheep loop tag, the SNAPP tag is lightweight yet durable and contains a waterproof and shock resistant transponder. This one piece tag locks together tightly without pinching the ear.

EiD TAGS ARE AVAILABLE IN ADDITIONAL COLOURS



SEE PAGE 68 FOR APPLICATORS

EiD BREEDING PAIRS

BUBBLEGUM

- A great 'tag for life' which is manufactured from an ultra-flexible plastic
- Concealing the transponder in the tag pin allows for excellent weight distribution
- The tip of the tag pin has a cross cut point which helps tag application by cutting through ear tissue as opposed to punching through
- Can be cut to create a two piece tag
- Lightweight yet robust applicator
- Management print options are available - see page 68
- **Applicator:** B49
- **Available colours (EiD)*:** Yellow, Blue, Green, Orange, Pink, Red and White
- **Available colours (visual):** Yellow, Blue, Green, Orange, Pink, Red, White, Brown, Grey, Purple and Spearmint



QWIK EID + QWIK Visual

- Manufactured from durable, high quality plastic designed for maximum retention
- Reliable EiD with superior read range transponder
- Anti-tamper hub
- Easy to read, fade resistant laser marking
- Tag point blunts on application for added safety
- High-performance cast aluminium, easy-load applicator
- Management print options are available - see page 68
- **Applicator:** 362
- **Available colours (EiD)*:** Yellow, Dark Blue, Dark Green, Orange, Red, White and Purple
- **Available colours (Visual):** Yellow, Light Blue, Dark Blue, Light Green, Dark Green, Orange, Pink, Red, White, Purple and Pink

* Depending on regional regulations

SHEEP TAGS

EiD BREEDING PAIRS

MINI EiD BUTTON

This new, lightweight EiD button tag enables flock managers to take full advantage of the management benefits of EiD; especially when tagging lambs at birth.

- Fits securely within the ear
- High quality laser marking
- Diameter: 24mm
- Choose from 3 visual match up options
- Management print options are available - see opposite
- **Applicator:** B08 (NB: please remove the black insert from the applicator before fitting the female EiD button)
- **Available colours*:**
Yellow, Blue, Green, Red and Purple



NEW TAG

USE THE SAME APPLICATOR AND PAIR WITH ONE OF THE FOLLOWING THREE VISUAL TAGS:

A. ULTRA BUTTON

Available colours: Yellow, Blue, Green, Orange, Pink, Red, White and Purple



ULTRA BUTTON



MINI EiD BUTTON

B. OVINA

Available colours: Yellow, Blue, Green, Orange, Pink, Red, White, Purple and Black



OVINA



MINI EiD BUTTON

C. ACE

Available colours: Yellow, Blue, Green, Orange, Pink, White, Purple and Black



ACE


MINI EiD BUTTON

* Depending on regional regulations

PEDIGREE/MANAGEMENT PRINT OPTIONS



QWIKvisual



BUBBLEGUM



OVINA



ACE



MINI EiD BUTTON



ULTRA BUTTON

**CONTACT
CUSTOMER
SERVICES TO
ORDER YOUR
FREE SAMPLES**

SHEEP AND GOATS

AVAILABLE EITHER PRINTED OR BLANK



BUBBLEGUM

MINIMUM ORDER 10

- **Applicator:** B49
- **Available colours:** Yellow, Blue, Green, Orange, Pink, Red, White, Brown, Grey, Purple and Spearmint



OVINA

NO MINIMUM ORDER

- **Applicator:** B08
- **Available colours:** Yellow, Blue, Green, Orange, Pink, Red, White, Purple and Black



MINI DALESMAN

MINIMUM ORDER 10

- **Applicator:** 09025
- **Available colours:** Yellow, Sky Blue, Lime Green, Orange, Pink, Red, White, Lilac and Grey



QWIKvisual

MINIMUM ORDER 10

- **Applicator:** 362
- **Available colours:** Yellow, Light Blue, Dark Blue, Light Green, Dark Green, Orange, Red, White, Purple and Pink



ADAMATIC

MINIMUM ORDER 10

- **Available Colours:** Yellow, Blue, Green, Orange, Pink, Red, White, Brown, Grey, Purple, Black, and Spearmint



D TAG

MINIMUM ORDER 10

- **Available Colours:** Blue, Green, Orange, Pink, Red, White, Purple and Black



ADAMATIC TAGGER

- **Product Code:** 66000103



D TAG TAGGER

- **Product Code:** B40

PIGS

AVAILABLE EITHER PRINTED OR BLANK



JUNIOR MEDIUM

- Large flag makes for easier visibility
- Flexible plastic
- Supplied blank or with high quality laser marking
- **Applicator:** B08
- **Available colours:** Yellow, Blue, Green, Pink, Red, Purple and White



OVINA

- Large flag makes for easier visibility
- Flexible plastic
- Supplied blank or with high quality laser marking
- **Applicator:** B08
- **Available colours:** Yellow, Blue, Green, Orange, Pink, Red, White, Purple and Black



STEEL PIG TAG

MINIMUM ORDER 10

- Made to withstand carcass heat treatment at slaughter
- Small steel tag, with curl lock
- Folded/crimped tag dimensions are 28mm x 6mm, with clearly etched ID
- **Applicator:** B45PIG
- Please note metal pig tags can take up to 10 working days to deliver



ULTRA BUTTON

- Diameter 28mm
- Good long term tag management infor can be printed on the back
- **Available Colours:** Yellow, Blue, Green Orange, Pink, Red and White

APPLICATORS

SPARE PINS AND INSERTS AVAILABLE
 Contact us for details

RapID SINGLE SHOT APPLICATOR
 PART CODE: **RAPIDEIDSINGLE**




RapID MULTISHOT APPLICATOR
 PART CODE: **RAPIDEIDMULTI**



BUBBLEGUM APPLICATOR
 PART CODE: **B49**



KLIK TAG APPLICATOR
 PART CODE: **ID02001**



SNAPP TAG APPLICATOR
 PART CODE: **09028**



QWIKID/QWIKvisual TAG APPLICATOR
 PART CODE: **362**



USE PIN TO
APPLY ULTRA
VISUAL TAGS



UNIVERSAL TISSUE
SAMPLE TAG APPLICATOR

PART CODE: B07MICRO




FIIP PIN APPLICATOR

PART CODE: B08



PIN
FLIPS
OUT



USE TO ADAPT YOUR
BLACK AND YELLOW TST
APPLICATOR FOR LIQUID
FILLED MICRO TUBES

PART CODE: 67001768



UNIVERSAL
APPLICATOR (RED)

PART CODE: B06



QWIKSILVER
APPLICATOR

PART CODE: 364



MINI DALESMAN
APPLICATOR

PART CODE: 09025



METAL PIG TAG
APPLICATOR

PART CODE: B45PIG



FIND OUT MORE - To browse the full range of our Livestock Identification products from the Allflex UK Group Portfolio please contact your local account manager or our customer services team.

A

Alligator Handling Systems	41
Antichew Liquid	38

B

Balling Gun	7
Bull Ring Applicator	13
Bull Rings	13

C

Calf Feeders	8
Calf Jacket	12
Calf Teats	8
Calving Aid	4
Calving Rope	4
Castration Ring Pliers	25
Castration Rings	24
Cattle Blow Dryer	57
Cattle/ Livestock Prod (Electric)	10
Cobalt Bolus	7
Colostrum (Calf & Lamb)	6 & 25
Cow Immobilizer (Kick Bar)	13
Cow Immobilizer Spring	13
Cow Immobilizer Tip	13
Crooks	26
Curry Combs	57

D

Dehorning Handles	10
Dehorning Shear Blades	11
Dehorning Wire	11
Digital Thermometer	34
Downacow Harness	12
Drenchers	9, 16, 17, 18 & 19

E

Emergency Feeders (Lamb)	25
Ewe Truss/ Prolapse Harness	23

F

Flutter Valve	19
Footcare & Cleansing Spray	20
Footrot Shears	20

G

Gambrel Restrainer	31
Gloves	5 & 30

H

Halters	14 & 57
Heavy Duty Hobbles	12
Heiniger Cattle Clippers	50, 51 & 52
Heiniger Clipper Blades	53 & 56
Heiniger Dog Clippers	51, 52 & 53
Heiniger Equine Clippers	50, 51 & 52
Heiniger Shearing Equipment	54 & 55
Heiniger Sheep Shear Sharpener	56
Hoof Cutter	20
Hoof Knives	21
Hooftrat Plus	20

I

Infrared Lamps	30
Infrared Reflector	30

K

Kamar Heatmount Detectors	37
---------------------------	----

L

Lamb Bar Feeder with Non Return Valve	29
Lamb Covers	31
Lamb Feeding Buckets & Bottles	28
Lamb Teats	29 & 30
Lambing Aid Set	23
Lambing Instrument	24
Lambing Rope	24
Livestock Identification	64 - 69

M

Marker Spray	31, 32 & 33
Mediclamp Clamp	34
Milk Sucking Preventers	13
Mobile Printer with Bluetooth	49

N

Needles	19 & 34
---------	---------

O

Obstetric Gel	4
Ovine Conditioner (Drench)	6

P

Pig Sorting Paddle	38
Pocket Knives	39
Prattley Animal Handling Systems	40

R

Raddle Powder	36
Ram Crayons	35 & 36
Ram Harnesses	35 & 36

S

S.W.A.T Auto Fly Control System	38
S.W.A.T Refill Can For Dispenser	38
Scissors	57
Sheep Marking Fluid	32
SI-RO-MARK	37
Speedy Feeder	9
Superstart Calf Feeding System	8
Superstart Stomach Tube	8
Suture Tape	33
Syringes	15 & 27

T

Tag Applicators	68 - 69
Tag Marker Pen	64
Tail Paint	36 & 37
Tail Tape	37
Tally Counter	55
Tar Plaster Spray	38
Trafalgar Shears	56
Tru - Test EID Readers	48
Tru -Test Suspension Cell	45
Tru -Test Electronic Weigh Indicators	43, 46, 47 & 49
Tru-Test Load Bars	44 - 45

U

Uteri Support	23
---------------	----

V

Vaccinators	15, 16, 18 & 27
-------------	-----------------

W

Wax Marker Stick	33
Whistles	38



Tel: **0845 6008081**

Fax: **0800 7836655**

Email: **tradesales@allflex.co.uk**

Web: **www.coxagri.com**
www.ritchey.co.uk

1 Greencroft Industrial Park
Stanley
County Durham
DH9 7YA

Cox and Ritchey are part of the Allflex UK Group Ltd

B17CR/TPB/09.17/SC
